

# m360 - PRODUCT CATALOGUE

---

06-2024



m360 SELECT SERIES

m360 MODIFICATIONS






MADE IN GERMANY.




m360

SOLUTIONS FOR YOU

# Content

	<b>m360 SELECT SERIES</b>	<b>4</b>
	CONSTANT CLIMATE CHAMBERS HPPCeco	<b>6</b>
	CLIMATE CHAMBERS ICH	<b>10</b>
	CLIMATE CHAMBERS ICHC	<b>13</b>
	UNIVERSAL OVENS U	<b>16</b>
	PASS-THROUGH OVENS UF TS	<b>19</b>
	VACUUM OVENS VOcool	<b>21</b>
	DISHWASHER UFdw	<b>23</b>
	PLASTICINE OVEN UFPIplus	<b>25</b>
	<b>m360 MODIFICATIONS</b>	<b>28</b>
	CONSTANT CLIMATE CHAMBERS HPPeco	<b>30</b>
	HUMIDITY CHAMBERS HCP	<b>33</b>
	CLIMATE CHAMBERS ICH	<b>34</b>
	ENVIRONMENTAL TEST CHAMBERS CTC / TTC	<b>36</b>
	UNIVERSAL OVENS U	<b>37</b>
	VACUUM OVENS VO	<b>40</b>
	INCUBATORS I	<b>41</b>
	CO <sub>2</sub> INCUBATORS ICO	<b>43</b>
	COMPRESSOR-COOLED INCUBATORS ICP	<b>44</b>
	PELTIER-COOLED INCUBATORS IPPeco	<b>46</b>
	<b>TECHNICAL INFORMATIONS MECHANICS</b>	<b>50</b>
	<b>TECHNICAL INFORMATIONS ELECTRONICS</b>	<b>72</b>
	<b>TECHNICAL INFORMATIONS ACCESSORY</b>	<b>92</b>



**m360**  
**SELECT SERIES**

CONSTANT CLIMATE CHAMBERS HPPCeco

CLIMATE CHAMBERS ICH

CLIMATE CHAMBERS ICHC

UNIVERSAL OVENS U

UNIVERSAL OVENS U

PASS-THROUGH OVENS UF TS

VACUUM OVENS VOcool


DISHWASHER UFdw

PLASTICINE OVEN UFPIplus



**CONSTANT CLIMATE CHAMBERS HPPCeco**

- Interior: stainless steel, mat. 1.4301 (ASTM 304), deep-drawn
- Housing: textured stainless steel 1.4016 (ASTM 430), rear zinc-plated steel 1.0037 (St37-2), intuitively operated TwinDISPLAY (TFT colour display) with touchscreen
- Door: Sizes 110 to 1060: fully insulated stainless steel doors with 2-point locking (compression door lock) and inner glass door, Sizes 1400 and 2200: stainless steel doors with glass sectors, fully heated inner glass panes integrated in the full-sight glass door with 2-point locking (compression door lock), sizes 750, 1060 and 1400 two leaves, size 2200 three leaves
- Connection: mains cable with plug (German type)
- Installation: Sizes 110 and 260: 4 adjustable appliance feet, Sizes 410 to 1060: on lockable castors, Sizes 1400 and 2200: on height-adjustable and lockable castors



Model sizes/Description			110	260	410	750	1060	1400	2200	
Stainless steel interior	Volume	l	108	256	384	749	1060	1360	2140	
	Width	(A) mm	560	640		1040		1250	1972	
	Height	(B) mm	480	800	1200		1450			
	Depth (less 32 mm for fan - Peltier)	(C) mm	400	500		600	850	750		
	Max. number of grids/shelves	number	5	9	14		28		42	
	Max. loading per grid/shelf	kg	20		30	60	30			
	Max. loading of chamber	kg	150	200				250	330	
	Max. loading per slide-in drip tray	kg	3	4		8		-		
	Max. loading per bottom drip tray	kg	3	4		8		-		
Textured stainless steel exterior	Width	(D) mm	745	824		1224		1435	2157	
	Height (sizes 410, 750, 1060, 1400 and 2200 with castors)	(E) mm	864	1183	1720		1913			
	Depth (without door handle, depth of handle +56 mm, +111 mm Peltier element)	(F) mm	555	655		755	1005	905		
Standard equipment	Standard works calibration certificate (measuring point chamber center)		+25 °C / 60 % rh, +40 °C / 75 % rh							
	Standard works calibration certificate (measuring point chamber center)	°C	5 %CO <sub>2</sub>							
	Stainless steel grid, electropolished	number	2					4	6	
	Water tank including connection hose		●							
	Door-open-recognition		●							

Model sizes/Description			110	260	410	750	1060	1400	2200	
Temperature	Working temperature range without light, without humidity	°C	0 (at least 20 below ambient temperature) up to +70							
	Working temperature range without light, with humidity	°C	+5 (at least 20 below ambient temperature) up to +70					+15 (at least 10 below ambient temperature) up to +60		
	Working temperature range with light, without or with humidity	°C	+15 to +40				-			
	Setting temperature range without light, with humidity	°C	+5 to +70					+15 to +60		
	Setting temperature range with light, with humidity	°C	+5 to +70				-			
	Setting temperature range with light, without humidity	°C	0 to +70				-			
	Setting temperature range without light, without humidity	°C	0 to +70							
	Setting accuracy	°C	0.1							
Humidity	Setting range humidity without light	% rh	10 to 90					10 to 80		
	Setting range humidity with light	% rh	10 to 90				-			
	Setting accuracy	% rh	0.5							
CO <sub>2</sub>	Setting range CO <sub>2</sub>	% CO <sub>2</sub>	0 to 20							
	Setting accuracy CO <sub>2</sub>	% CO <sub>2</sub>	0.1							
	Variation in time CO <sub>2</sub>	% CO <sub>2</sub>	± 0,2							
Further data	Electrical load at 230 V, 50/60 Hz	W	420	700		1400		2000		
	Electrical load at 115 V, 50/60 Hz	W	420	700		1400		-		
	Peltier elements in the rear	number	1	2		4		6		
Packing data	Net weight	kg	64	101	142	213	259	347	474	
	Gross weight (packed in carton)	kg	87	140	192	279	339	525	718	
	Width	mm	830	930		1330	1370	1560	2300	
	Height	mm	1050	1380	1930	1910	1970	2200		
	Depth	mm	800	930		1050	1300	1190		

**Order No. Constant Climate Chambers**      HPP110Ceco | HPP260Ceco | HPP410Ceco | HPP750Ceco | HPP1060Ceco | HPP1400Ceco | HPP2200Ceco

Options		110	260	410	750	1060	1400	2200	
Temperature restriction; Temperature: +60 °C		A8						-	
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber - standard with 1060) - includes replacement of standard grids by reinforced grids		-				K1		-	
Light module cold white 6,500 K: LED light strips arranged on the side walls of the interior, 10 for model 110, 14 for model 260/410/750, programme controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature and humidity		T7				-			
Light module cold white 6,500 K + warm white 2,700 K: LED light strips - 10 stripes for model 110, 14 for model 260/410/750 - (5 resp. 7 alternating cold white light strips and 5 resp. 7 warm white light strips) on the side walls of the interior, programme-controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature and humidity		T8				-			
Light module warm white 2,700 K: light strips arranged on the side walls of the interior, 10 strips for model 110, 14 for model 260/410/750, programme-controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature and humidity		T9				-			
Interior socket (can only be ordered with limited temperature range - max. +60 °C), ampacity 230 V / 2,2A / 500W, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68 (requires option A8 against surcharge) (max. one per appliance)		R3				-			
Entry port with 23 mm clear diameter, can be closed, standard position (F0 and F2 not for model size 260 with light module; F0 - F3 not for model size 110 with light module)	left centre/centre	F0				-			
	left centre/top	F1				-			
	right centre/centre	F2				-			
	right centre/top	F3				-			



Options		110	260	410	750	1060	1400	2200
Entry port with 23 mm clear diameter, can be closed, individual position (please state position)	left				F4			
	right				F5			
	rear				F6			-
Entry port (silicone), 40 mm clear diameter, moisture tight, can be closed by silicone stopper, at the back (please state location)					F7			-
Entry port (silicone), 40 mm clear diameter, moisture tight, can be closed by silicone stopper, standard position rear					D8			
4 - 20 mA current loop interface	Temperature controller actual value (-10 to +80 °C = 4 - 20 mA)				V3			
	CO <sub>2</sub> controller, actual value (0 to 25 % CO <sub>2</sub> = 4 - 20 mA)				V9			
	Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring - price per sensor				V6			
	Humidity controller, actual value (0 to 100 % rh = 4 - 20 mA)				V7			
Compressed air dehumidification (efficient dehumidification of the interior by means of compressed air). Standard works calibration certificate (measuring point chamber centre) at +10 °C with 10 % rh					C9			-
Door with lock and key (safety lock)	One lock				B6			-
	Two locks (one each door)				-	B62		-
	Three locks (one each door)				-			B63
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)					H6			
Process-dependent programmable door lock	One				D4			-
	Two (one each door)				-	D42		-
	Three (one each door)				-			D43
Door-open-recognition	One				V5			-
	Two (one each door)				-	V52		-
	Three (one each door)				-			V53
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)					C3			
MobileALERT for 2 alarm notifications; notification by SMS. temperature and humidity alarm (requires option H6)					C4			
Castor frame (2-part), height 140 mm		R9						

Accessories	110	260	410	750	1060	1400	2200
Stainless steel grid, electropolished	E20165	E28891		E20182	B41251	B38955	
Reinforced stainless steel grid, electropolished, max. loading 60 kg; size 750 with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber	E29767	E29766		B32190	B32550	-	
Perforated stainless steel shelf	B00325	B29725		B00328			-
Reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber				B32191	B32549	-	
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (may affect the temperature distribution, requires option K1)				B32763			-
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1 for sizes 450/750)	E02073	B29726		E02075	B32599	-	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, requires option K1)				B34055			-
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1 for sizes 450/750)	B04359	B29722		B04362	B29769	-	
Holder for water tank (sizes 110 - 750: 2.5 litres, sizes 1060/1400/2200: 10 litres) for mounting on the rear of the appliance. Standard equipment for sizes 410, 750, 1060, 1400 and 2200	B49092					-	
Central water supply with filter cartridges for connection to the domestic water supply. Product information on demand				ZWVR6			
Central water supply without filter cartridges for connection to the domestic water supply (only for demineralised water with a conductivity of 5 to 10 µS/cm and a pH value between 5 and 7). Product information on demand				ZWVR7			




Accessories	110	260	410	750	1060	1400	2200
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29734	B29738	B42116	B29742	-	-	-
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29735	B29739	B42117	B29743	-	-	-
Stacking set (4 pcs) for stacking of appliances of same size	B29744	-	-	-	-	-	-
Set of height adjustable feet (4 pcs)	B29768	-	-	-	-	-	-
Subframe, adjustable in height (height 500 mm)	B29749	B29751	-	-	-	-	-
Subframe, on castors (height 560 mm)	B29750	-	-	-	-	-	-
Subframe, adjustable in height, height 130 mm	B33661	B33664	-	-	-	-	-
External measuring instrument with sensors for daylight and UV-light. Product information on demand	B04713				-	-	-
External measuring instrument with additional measuring head for temperature and humidity measurement. Product information on demand	B04714				-	-	-

Accessories   Software and Connectivity	110	260	410	750	1060	1400	2200
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Respective IQ/OQ documents available in German and English language (without surcharge)				FDAQ1			
Integration of additional units (up to max. 31 units) into an already existent FDA-software licence				FDAQ2			
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number				B33170			
Ethernet connection cable 5 m for computer interface				E06189			
USB-Ethernet adapter				E06192			

**CLIMATE CHAMBERS ICH**

- Interior: stainless steel, mat. 1.4301 (ASTM 304), deep-drawn
- Housing: textured stainless steel 1.4016 (ASTM 430), rear zinc-plated steel 1.0037 (St37-2), intuitively operated TwinDISPLAY (TFT colour display) with touchscreen
- Door: fully insulated stainless steel doors with 2-point locking (compression door lock) and inner glass door, from size 450 two leaves
- Connection: mains cable with plug (German type)
- Installation: on lockable castors
- Interfaces:
 

USB	Ethernet
	LAN



Model sizes/Description			450	450L
Stainless steel interior	Volume	l		449
	Width	(A) mm		1040
	Height	(B) mm		720
	Depth (less 33 mm for fan)	(C) mm		600
	Max. number of grids/shelves	number		14
	Max. loading per grid/shelf	kg		30
	Max. loading of chamber	kg		200
	Max. loading per slide-in drip tray	kg		8
Max. loading per bottom drip tray	kg		8	
Textured stainless steel exterior	Width	(D) mm		1224
	Height (with castors)	(E) mm		1224
	Depth (without door handle, depth of handle +56 mm)	(F) mm		785
Standard equipment	Standard works calibration certificate (measuring point chamber center)	°C	+10 °C , +37 °C , +30 °C / 60 % rh	
	Stainless steel grid, electropolished	number		2
	Entry port (silicone), 40 mm clear diameter, moisture tight, can be closed by a silicone stopper, standard position at the back			●
	Water tank including connection hose			●
	Door-open-recognition		-	●
Temperature	Working temperature range without humidity ICH	°C	-10 to +60	
	Working temperature range ICH /ICH L with humidity and/or light	°C	+10 to +60	
	Working temperature range ICH L without humidity	°C	0 to +60	
	Setting temperature range ICH	°C	-10 to +60	
	Setting temperature range ICH L	°C	0 to +60	
	Setting accuracy	°C	0.1	
Humidity	Setting range humidity	% rh	10 to 80	
	Setting accuracy	% rh	0.5	
Light	Illumination unit (only ICH L) acc. ICH Q1B, option 2; separately switchable via controller, one box; number of fluorescent lights with cold white light (size 110: 3, size 260/750: 4), number of fluorescent lights with UV lamps (all sizes: 2)		-	●
Further data	Electrical load at 230/115 V, 50/60 Hz ICH	W	1350	1550

Model sizes/Description			450	450L
Packing data	Net weight	kg	220	
	Gross weight (packed in carton)	kg	285	
	Width	mm	1330	
	Height	mm	2150	
	Depth	mm	1700	
<b>Order No. Climate Chambers</b> ICH = Climate chamber ICH L = Climate chamber with light			ICH450	ICH450L



Options			450	450L	
Voltage 115 V, 50/60 Hz			X2		
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) - includes replacement of standard grids by reinforced grids (ICH)			K1	-	
Alternative light boxes (replace the standard lighting; have to be ordered together with the chamber); number of fluorescent lamps: size 110: 5, sizes 260/750: 6, with cold white light according to ICH guideline Q1B, option 2 (only ICH L)			one box	T81	
Alternative light boxes (replace the standard lighting; have to be ordered together with the chamber); number of fluorescent UV lamps: size 110: 5, sizes 260/750: 6, with UV-A radiation according to ICH guideline Q1B, option 2 (only ICH L)			one box	T01	
Interior socket, ampacity 230V / 2,2A / 500W, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68 (not for ICH110L) (max. one per appliance)			R3		
Entry port with 23 mm clear diameter, can be closed, standard position (F1 and F3 not for models ICH L)			left centre/centre	F0	-
			left centre/top	F1	-
			right centre/top	F3	-
Entry port (silicone), 40 mm clear diameter, moisture tight, can be closed by silicone stopper, at the back (please state location). Not for models ICH L			F7	-	
4 - 20 mA current loop interface			Humidity controller, actual value (0 to 100 % rh = 4 - 20 mA)		V7
			Temperature controller actual value (-20 to +70 °C = 4 - 20 mA)		V3
			Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 3) - price per sensor (-20 to +70 °C = 4 - 20 mA)		V6
Fan speed monitoring with switching off the heating and with alarm in case of failure			V4		
Compressed air dehumidification (efficient dehumidification of the interior by means of compressed air). Standard works calibration certificate (measuring point chamber centre) at +10 °C with 10 % rh			C9		
Door with lock and key (safety lock)			B6		
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)			H6		
Door-open-recognition, shuts down humidity, light (standard with ICH L)			V5	-	
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)			C3		


Accessories			450	450L
Stainless steel grid, electropolished			E20182	
Reinforced stainless steel grid, electropolished, max. loading 60 kg; size 750 with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber			B32190	-
Perforated stainless steel shelf			B00328	
Reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber			B32191	-
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (may affect the temperature distribution, requires option K1)			B32763	
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1 for sizes 450/750)			E02075	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1 for sizes 450/750)			B04362	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, requires option K1)			B34055	
Central water supply with filter cartridges for connection to the domestic water supply. Product information on demand			ZWVR6	
Central water supply without filter cartridges for connection to the domestic water supply (only for demineralised water with a conductivity of 5 to 10 µS/cm and a pH value between 5 and 7). Product information on demand			ZWVR7	
External measuring instrument with sensors for daylight and UV-light. Product information on demand (only ICH L)			-	B04713
External measuring instrument with additional measuring head for temperature and humidity measurement. Product information on demand			B04714	



Accessories   Software and Connectivity	450	450L
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1	
Integration of additional units (up to max. 31 units) into an already existent FDA-software licence	FDAQ2	
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number	B33170	
Ethernet connection cable 5 m for computer interface	E06189	
USB-Ethernet adapter	E06192	

CLIMATE CHAMBERS ICHC

- Interior: stainless steel, mat. 1.4301 (ASTM 304), deep-drawn
- Housing: textured stainless steel 1.4016 (ASTM 430), rear zinc-plated steel 1.0037 (Sf37-2), intuitively operated TwinDISPLAY (TFT colour display) with touchscreen
- Door: fully insulated stainless steel doors with 2-point locking (compression door lock) and inner glass door, size 750: two leaves
- Connection: mains cable with plug (German type)
- Installation: on lockable castors
- Interfaces:  



Model sizes/Description			110	260	750
Stainless steel interior	Volume	l	108	256	749
	Width	(A) mm	560	640	1040
	Height	(B) mm	480	800	1200
	Depth (less 33 mm for fan)	(C) mm	400	500	600
	Max. number of grids/shelves	number	5	9	14
	Max. loading per grid/shelf	kg	20		30
	Max. loading of chamber	kg	150	200	
	Max. loading per slide-in drip tray	kg	3	4	8
	Max. loading per bottom drip tray	kg	3	4	8
Textured stainless steel exterior	Width	(D) mm	745	824	1224
	Height (with castors)	(E) mm	1233	1552	1950
	Depth (without door handle, depth of handle +56 mm)	(F) mm	585	685	785
Standard equipment	Standard works calibration certificate (measuring point chamber center)		+10 °C , +37 °C , +30 °C / 60 % rh		
	Standard works calibration certificate (measuring point chamber center)	°C	5 %CO <sub>2</sub>		
	Stainless steel grid, electropolished	number	2		
	Entry port (silicone), 40 mm clear diameter, moisture tight, can be closed by a silicone stopper, standard position at the back		●		
	Water tank including connection hose		●		
	CO <sub>2</sub> connection set: hose with coupling and clamp		●		
	Door-open-recognition		●		
Temperature	Working temperature range without humidity	°C	0 to +60		
	Working temperature range with humidity	°C	+10 to +60		
	Setting temperature range	°C	0 to +60		
	Setting accuracy	°C	0.1		
Humidity	Setting range humidity	% rh	10 to 80		
	Setting accuracy	% rh	0.5		
CO <sub>2</sub>	Setting range CO <sub>2</sub>	%CO <sub>2</sub>	0 to 20		
	Setting accuracy CO <sub>2</sub>	%CO <sub>2</sub>	0,1		
	Variation in time CO <sub>2</sub>	%CO <sub>2</sub>	± 0,2		
Further data	Electrical load at 230/115 V, 50/60 Hz ICH	W	1350		



Model sizes/Description			110	260	750
Packing data	Net weight	kg	114	160	249
	Gross weight (packed in carton)	kg	142	217	319
	Width	mm	880	930	1330
	Height	mm	1410	1760	2150
	Depth	mm	810	930	1050
<b>Order No. Climate Chambers</b> ICH = Climate chamber ICH L = Climate chamber with light			ICH110C	ICH260C	ICH750C

Options		110	260	750
Voltage 115 V, 50/60 Hz		X2		
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) - includes replacement of standard grids by reinforced grids (ICH)		-		K1
Entry port with 23 mm clear diameter, can be closed, standard position (F1 and F3 not for models ICH L)	left centre/centre	F0		
	left centre/top	F1		
	right centre/top	F3		
Entry port (silicone), 40 mm clear diameter, moisture tight, can be closed by silicone stopper, at the back (please state location)		-	F7	
4 - 20 mA current loop interface	Temperature controller actual value (-20 to +70 °C = 4 - 20 mA)	V3		
	Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 3) - price per sensor (-20 to +70 °C = 4 - 20 mA)	V6		
	Humidity controller, actual value (0 to 100 % rh = 4 - 20 mA)	V7		
	CO <sub>2</sub> controller, actual value (0 to 25 % CO <sub>2</sub> = 4 - 20 mA)	V9		
Fan speed monitoring with switching off the heating and with alarm in case of failure		V4		
Compressed air dehumidification (efficient dehumidification of the interior by means of compressed air). Standard works calibration certificate (measuring point chamber centre) at +10 °C with 10 % rh		C9		
Door with lock and key (safety lock)		B6		
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)		H6		
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)		C3		

Accessories		110	260	750
Stainless steel grid, electropolished		E20165	E28891	E20182
Reinforced stainless steel grid, electropolished, max. loading 60 kg; size 750 with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber		E29767	E29766	B32190
Perforated stainless steel shelf		B00325	B29725	B00328
Reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber		-		B32191
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (may affect the temperature distribution, requires option K1)		-		B32763
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1 for sizes 450/750)		E02073	B29726	E02075
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, requires option K1)		-		B34055
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1 for sizes 450/750)		B04359	B29722	B04362
Holder for water tank (2.5 litres) for mounting on the rear of the appliance. Standard equipment for size 750		E32172		-
Central water supply with filter cartridges for connection to the domestic water supply. Product information on demand		ZWVR6		
Central water supply without filter cartridges for connection to the domestic water supply (only for demineralised water with a conductivity of 5 to 10 µS/cm and a pH value between 5 and 7). Product information on demand		ZWVR7		
External measuring instrument with additional measuring head for temperature and humidity measurement. Product information on demand		B04714		


Accessories   Software and Connectivity		110	260	750
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Respective IQ/OQ documents available in German and English language (without surcharge)		FDAQ1		
Integration of additional units (up to max. 31 units) into an already existent FDA-software licence		FDAQ2		
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number		B33170		

Accessories   Software and Connectivity	110	260	750
Ethernet connection cable 5 m for computer interface		E06189	
USB-Ethernet adapter		E06192	

**UNIVERSAL OVENS U**

- Interior: stainless steel, material 1.4301 (ASTM 304), with all-round deep-drawn ribs to integrate the large-area heating with ceramic-metal sheath
- Housing: textured stainless steel 1.4016 (ASTM 430), rear zinc-plated steel 1.0037 (St37-2), intuitively operated SingleDISPLAY or TwinDISPLAY (TFT colour display) with touchscreen
- Door: fully insulated stainless steel doors with 2-point locking (compression door lock), from size 450 two leaves
- Connection: mains cable with plug (German type) CEE plug for 400 V
- Installation: Size 250: 4 adjustable appliance feet, Sizes 400 to 2500: on lockable castors
- Interfaces:

USB: only TwinDISPLAY



Model sizes/Description			250	400	600	1550	2500
Stainless steel interior	Volume	l	256	384	600	1550	2554
	Width	(A) mm	800	640		1040	1680
	Height	(B) mm	640	1200	1680	1520	
	Depth (less 39 mm for fan)	(C) mm	500		1000		
	Max. number of grids/shelves	number	7	14	20	18	
	Max. loading per grid/shelf	kg	20	30	20	30	60
	Max. loading of chamber	kg	300				
	Max. loading per slide-in drip tray	kg	4	8			
	Max. loading per bottom drip tray	kg	4	8			
Textured stainless steel exterior	Width	(D) mm	984	824		1224	1861
	Height (sizes 450, 750, 1060 with castors)	(E) mm	1023	1720	2200	1982	1983
	Depth (without door handle, depth of handle +56 mm)	(F) mm	684		1186		
Standard equipment	Standard works calibration certificate (measuring point chamber center)	°C	+ 160				
	Stainless steel grid, electropolished	number	2				
Temperature	Working temperature range	°C	at least 5 (UNplus) 10 (UFplus) above ambient temperature up to +300				
	Setting temperature range	°C	+20 to +300				
	Setting accuracy	°C	up to 99.9: 0.1 / from 100: 0.5				
Further data	Electrical load at 230 V, 50/60 Hz	W	3350		-		
	Electrical load at 400 V and 3 x 230V w/o neutral, 50/60 Hz	W	-	6500	7000	8000	
Packing data	Net weight	kg	110	161	135	217	675
	Gross weight (packed in carton)	kg	161	227	296	288	1020
	Width	mm	1090	930		1336	1980
	Height	mm	1360	1930	2320	2264	2260
	Depth	mm	930		1430		
<b>Order No. Universal Ovens</b>			UN250plus	UN400plus	-		
U = Universal Oven							
N = Natural convection							
F = Forced air circulation							
plus = Model with TwinDISPLAY			UF250plus	UF400plus	UF600plus	UF1550plus	UF2500plus



Options	250	400	600	1550	2500
Voltage 115 V, 50/60 Hz	X2		-		
Temperature restriction; Temperatures: +60, +70, +80, +95, +100, +120, +160, +180, +200 or +250 °C (Please indicate upon ordering)	A8				
Full-sight glass door (4-layer insulating glass) Temperature-range up to max. +250 °C (requires option A8 against surcharge)	B0		-		
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) - includes replacement of standard grids by reinforced grids	K1		-	K1	
Interior lighting for observing the load	R0				
Interior socket (can only be ordered with limited temperature range - max. +60 °C), ampacity 230 V / 2,2A / 500W, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68 (requires option A8 against surcharge) (max. one per appliance)	R3				
Interior nearly gastight - technical clarification required for combination with other options	K2			-	
Interior nearly gastight with possibility for gas inlet/outlet through 2 tubes with ball valves - technical clarification required for combination with other options	K3			-	
Entry port with 23 mm clear diameter, can be closed, standard position	left centre/centre	F0			
	left centre/top	F1			
	right centre/centre	F2			
	right centre/top	F3			
Entry port with 23 mm clear diameter, can be closed, individual position (please state position)	left	F4			
	right	F5			
	rear	F6			
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please state location)	D6				
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please state location)	F7				
Entry port with 57 mm clear diameter, can be closed by flap, individual positions in the back wall (please state position, not possible for devices with forced air circulation of sizes 30-75)	F8				
Entry port with 100 mm clear diameter, can be closed by flap, individual positions in the back wall (please state position, not possible for devices with forced air circulation of size 75)	F9				
4 - 20 mA current loop interface (0 to +310 °C = Temperature controller, actual value 4 - 20 mA) Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 1 SingleDISPLAY, max. 3 TwinDISPLAY) - price per sensor	V3				
	V6				
Fan speed monitoring with switching off the heating and with alarm in case of failure (for UFplus)	V4		-		
Door with lock and key (safety lock)	B6				
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)	H6				
Potential-free contact (24 V/2 A) with socket, for signal generation, controlled by programme segment, for free-selectable functions to be activated (e.g. activation of audible and visual signals, exhaust motors, fans, stirrers, etc.). Only for units with TwinDISPLAY	2 contacts		H72		
Process-dependent programmable door lock (only for units with TwinDISPLAY)	D4				
Door-open-recognition (only for units with TwinDISPLAY)	V5				
Flexible Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 1 additional sensor is possible). The measured temperature will be indicated on the display, recorded in the integral data store, and can be documented via the AtmoCONTROL software	H8				
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)	C3				
Castor frame (2-part), height 140 mm	R9	-			
Fresh-air filter (filtration efficiency 80 %) mounted at the device bottom (for devices with forced air circulation, sizes 30 - 260 require castor frame R9 or subframe)	R8			-	



Accessories	250	400	600	1550	2500
Reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber	-			B50539	B58583
Reinforced stainless steel grid, electropolished, max. loading 60 kg; from size 450 with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber	-	E29766		B58238	B57167
Stainless steel grid, electropolished	E20167	E28891		-	
Perforated stainless steel shelf	B00326	B29725		-	
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution)	-	B29726		-	




Accessories	250	400	600	1550	2500
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution)	-	B29722		-	
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), straight, for exhaust air ducting (if necessary for connection by hose)	B29718				
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose)	B29719				
Set of height adjustable feet (4 pcs)	B29768	-			

Accessories   Software and Connectivity	250	400	600	1550	2500
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1				
Integration of one additional unit (up to max. 31 units) into an already existent FDA-software licence (only for units with TwinDISPLAY)	FDAQ2				
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)	B33170				
Ethernet connection cable 5 m for computer interface	E06189				
USB-Ethernet adapter	E06192				

**PASS-THROUGH OVENS UF TS**

- Interior: stainless steel, material 1.4301 (ASTM 304), with all-round deep-drawn ribs to integrate the large-area heating with ceramic-metal sheath
- Housing: textured stainless steel 1.4016 (ASTM 430), intuitively operated TwinDISPLAY (TFT colour display) with touchscreen
- Door: fully insulated stainless steel doors with 2-point locking (compression door lock)
- Connection: mains cable with plug (German type) CEE plug for 400 V
- Installation: 4 adjustable appliance feet
- Interfaces:  



Model sizes/Description				1250
Stainless steel interior	Volume	l		1248
	Width	(A) mm		1040
	Height	(B) mm		1200
	Depth	(C) mm		1000
	Max. number of grids/shelves	number		14
	Max. loading per grid/shelf	kg		30
	Max. loading of chamber	kg		300
	Max. loading per slide-in drip tray	kg		8
Max. loading per bottom drip tray	kg		8	
Textured stainless steel exterior	Width	(D) mm		1224
	Height	(E) mm		1714
	Depth (without door handle, depth of handle 2x +56 mm)	(F) mm		1181
Standard equipment	Standard works calibration certificate (measuring point chamber center)	°C		+ 160
	Stainless steel grid, electropolished	number		2
Temperature	Working temperature range	°C		at least 10 above ambient temperature up to +250
	Setting temperature range	°C		+20 to +250
	Setting accuracy	°C		up to 99.9: 0.1 / from 100: 0.5
Further data	Electrical load at 400 V and 3 x 230V w/o neutral, 50/60 Hz	W		7000
Packing data	Net weight	kg		372
	Gross weight (packed in carton)	kg		600
	Width	mm		1330
	Height	mm		1920
	Depth	mm		1050
<b>Order No. Pass-Through Ovens</b>				<b>UF1250TS</b>






Options	1250	
Temperature restriction; Temperatures: +60, +70, +80, +95, +100, +120, +160, +180 or +200 °C (Please indicate upon ordering)	A8	
Full-sight glass door (4 layer insulating glass) - extra cost per side - Temperature-range up to max. +250 °C	B0	
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) - includes replacement of standard grids by reinforced grids	K1	
Entry port with 23 mm clear diameter, can be closed, standard position	left centre/centre	F0
	left centre/top	F1
	right centre/centre	F2
	right centre/top	F3
	Entry port with 23 mm clear diameter, can be closed, individual position (please state position)	left
right		F5
Locking mechanism to prevent simultaneous opening of doors for contamination protection in case of wall installation (requires option V5)	D5	
4 - 20 mA current loop interface (0 to +260 °C = 4 - 20 mA)	Temperature controller, actual value	V3
	Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 3) - price per sensor	V6
Door with lock (safety lock); per side	B6	
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)	H6	
Potential-free contact (24 V/2 A) with socket, for signal generation, controlled by programme segment, for free-selectable functions to be activated (e.g. activation of audible and visual signals, exhaust motors, fans, stirrers, etc.). Only for units with TwinDISPLAY	2 contacts H72	
Process-dependent electromagnetic door lock (both sides)	D4	
Door-open-recognition; per side	V5	
Flexible Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 1 additional sensor is possible). The measured temperature will be indicated on the display, recorded in the integral data store, and can be documented via the AtmoCONTROL software	H8	
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)	C3	

Accessories	1250
Stainless steel grid, electropolished	B58903
Perforated stainless steel shelf	B58948
Flush-fit unit set (stainless steel frame covering gap between oven and wall opening), without air slots - technical clarification required	B33207

Accessories   Software and Connectivity	1250
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1
Integration of additional units (up to max. 31 units) into an already existent FDA-software licence	FDAQ2
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number	B33170
Ethernet connection cable 5 m for computer interface	E06189
USB-Ethernet adapter	E06192

VACUUM OVENS VOcool

- Interior: stainless steel interior, material 1.4404 (ASTM 316 L), hermetically welded, with removable mountings at the sides for cleaning, including thermoshelf guide bars, as well as mounting on top to avoid turbulences
- Housing: textured stainless steel 1.4016 (ASTM 430), rear zinc-plated steel 1.0037 (Si37-2), intuitively operated TwinDISPLAY (TFT colour display) with touchscreen
- Door: full-sight glass door, inside spring-loaded, 1.5mm thick glazed panel in safety glass, outside with anti-splitter screen
- Connection: mains cable with plug (German type)
- Installation: 4 adjustable appliance feet
- Interfaces:  



Model sizes/Description			29	49	101
Stainless steel interior	Volume	l	29	49	101
	Width	(A) mm	385		545
	Height	(B) mm	305	385	465
	Depth	(C) mm	250	330	400
	Max. number of thermoshelves	number	1		
	Max. loading of chamber	kg	40	60	
	Max. loading per thermoshelf	kg	20		
Textured stainless steel exterior	Width	(D) mm	550		710
	Height	(E) mm	607	687	767
	Depth (without door handle, depth of handle +38 mm)	(F) mm	400	480	550
Standard equipment	Standard works calibration certificate		+ 80 °C / 20 mbar		
	Thermoshelves - aluminium eloxadised, mat. 3.3547 (ASTM B209) - with integrated large-area heating including local temperature sensing (Pt100, 4-wire-circuit); individual overtemp. protection for each shelf. Further data see stainless steel number inner working chamber	number	1		
Temperature	Working temperature range	°C	+5 to +90		
	Setting temperature range	°C	+5 to +90		
	Setting accuracy	°C	0.1		
	Temperature variation in time (aluminium thermoshelf)	K	≤ ± 0,3		
	Temperature uniformity (surface) at +80 °C/20 mbar (aluminium thermoshelf)	K	≤ ± 0,5		
Pressure (vacuum)	Adjustment range pressure	mbar	5 to 1100		
	Permitted final vacuum	mbar	0.01		
	Maximum leakage rate	bar/h	0.01		
Further data	Electrical load at 230 V, 50/60 Hz	W	500	600	750
Packing data	Net weight	kg	66	95	130
	Gross weight (packed in carton)	kg	87	116	160
	Packed dimensions (Width, Height, Depth)		660/870/590		830/1050/800
	Net weight pump module without/with pump	kg	25/41	30/46	41/57
	Gross weight pump module without/with pump (packed in carton)	kg	46/62	51/67	66/82
	Packed dimensions pump module (Width, Height, Depth)	mm	660/870/590		830/1050/800
<b>Order No. Vacuum Ovens</b>			VO29cool	VO49cool	VO101cool




Options	29	49	101
Premium module: comprises the inert gas inlet, extra connectors for thermoshelves, 1 (size 29), 2 (sizes 49/101)	T5		
4 - 20 mA current loop interface (requires option T5)	Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring - price per sensor	V6	
		Temperature controller, actual value (0 to +100 °C = 4 - 20 mA)	
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)	H6		
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)	C3		
MobileALERT for 2 alarm notifications; notification by SMS. temperature and vacuum alarm (requires option T5 and H6)	-	C4	

Accessories	29	49	101
Subframe, tubular steel, black enamelled (for stacking unit consisting of vacuum oven and pump module, total height: 1650 mm, see „further data“ and sketch of oven dimensions)	E02030	E02031	E02037
Signal cable (3 m) for control of rotation speed and optimising pump performance by demand-controlled activation of purge of Memmert pump (not required with pump module)	B39410		
Noise-insulated vacuum pump module without pump (exterior dimensions and -material No. s. vacuum oven) with anti-vibration metal plate at the bottom to accommodate the vacuum pump, incl. full-sight glass door. Socket, signal cable and connecting hose to the vacuum oven	PM29	PM49	PM101
Noise-insulated vacuum pump module, as above, however with built-in pump, 230 V, 50/60 Hz	PMP29	PMP49	PMP101
Vacuum connecting hose (3 m) from oven to Memmert pump incl. optimised connection accessories (partially stainless steel), (not required with pump module)	B04026		
Vacuum pump with chemically resistant 4x diaphragm, pump capacity at atm. pressures: approx. 50 NI./min = 3,0 m³/h and autom. purge control, 230 V, 50 Hz. Max. guarantee period 2 years (requires accessories B39410 and B04026)	E07509		


Accessories   Software and Connectivity	29	49	101
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1		
Integration of additional units (up to max. 31 units) into an already existent FDA-software licence	FDAQ2		
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number	B33170		
Ethernet connection cable 5 m for computer interface	E06189		
USB-Ethernet adapter	E06192		

## DISHWASHER UFdw

- Interior: stainless steel, material 1.4301 (ASTM 304), with all-round deep-drawn ribs to integrate the large-area heating with ceramic-metal sheath
- Housing: textured stainless steel 1.4016 (ASTM 430), rear zinc-plated steel 1.0037 (S137-2), intuitively operated SingleDISPLAY or TwinDISPLAY (TFT colour display) with touchscreen
- Door: fully insulated stainless steel doors with 2-point locking (compression door lock), size 750: two leaves
- Connection: mains cable with plug (German type) CEE plug for 400 V
- Installation: on lockable castors
- Interfaces:



USB: only TwinDISPLAY



Model sizes/Description			750
Stainless steel interior	Volume	l	749
	Width	(A) mm	1040
	Height	(B) mm	1200
	Depth (less 39 mm for fan)	(C) mm	600
	Max. number of grids/shelves	number	14
	Max. loading per grid/shelf	kg	30
	Max. loading of chamber	kg	300
	Max. loading per slide-in drip tray	kg	8
	Max. loading per bottom drip tray	kg	8
Textured stainless steel exterior	Width	(D) mm	1224
	Height (sizes 450, 750, 1060 with castors)	(E) mm	1720
	Depth (without door handle, depth of handle +56 mm)	(F) mm	784
Standard equipment	Standard works calibration certificate (measuring point chamber center)	°C	+80
	CO <sub>2</sub> connection set: hose with coupling and clamp		●
Temperature	Working temperature range	°C	at least 10 above ambient temperature up to +100
	Setting temperature range	°C	+20 to +100
	Setting accuracy	°C	up to 99.9: 0.1 / from 100: 0.5
Further data	Electrical load at 400 V and 3 x 230V w/o neutral, 50/60 Hz	W	4800
Packing data	Net weight	kg	217
	Gross weight (packed in carton)	kg	288
	Width	mm	1330
	Height	mm	1910
	Depth	mm	1050

**Order No. Universal Ovens**

U = Universal Oven

N = Natural convection

F = Forced air circulation

plus = Model with TwinDISPLAY


UF750dw




Options		750
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)		H6
Potential-free contact (24 V/2 A) with socket, for signal generation, controlled by programme segment, for free-selectable functions to be activated (e.g. activation of audible and visual signals, exhaust motors, fans, stirrers, etc.). Only for units with TwinDISPLAY	2 contacts	H72
Process-dependent programmable door lock		D4
Door-open-recognition		V5
Accessories		750
Additional loading set consisting of 8 grids as well as 6 place settings and 3 grid baskets for cups/glasses, made of stainless steel, material 1.4301 (ASTM 304)		B51799
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), straight, for exhaust air ducting (if necessary for connection by hose)		B29718
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose)		B29719
Accessories   Software and Connectivity		750
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)		FDAQ1
Integration of one additional unit (up to max. 31 units) into an already existent FDA-software licence (only for units with TwinDISPLAY)		FDAQ2
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)		B33170
Ethernet connection cable 5 m for computer interface		E06189
USB-Ethernet adapter		E06192



**PLASTICINE OVEN UFPIplus**

- Interior:** stainless steel, material 1.4301 (ASTM 304), with all-round deep-drawn ribs to integrate the large-area heating with ceramic-metal sheath
- Housing:** textured stainless steel 1.4016 (ASTM 430), rear zinc-plated steel 1.0037 (S137-2), intuitively operated SingleDISPLAY or TwinDISPLAY (TFT colour display) with touchscreen
- Door:** fully insulated stainless steel doors with 2-point locking (compression door lock), from size 450 two leaves
- Connection:** mains cable with plug (German type) CEE plug for 400 V
- Installation:** Size 260: 4 adjustable appliance feet, Sizes 450 to 1060: on lockable castors
- Interfaces:**


USB: only TwinDISPLAY



Model sizes/Description			260	450	750	1060
Stainless steel interior	Volume	l	256	449	749	1060
	Width	(A) mm	640	1040		
	Height	(B) mm	800	720	1200	
	Depth (less 39 mm for fan)	(C) mm	500	600		850
	Max. number of grids/shelves	number	9	8	14	
	Max. loading per grid/shelf	kg	20	30		60
	Max. loading of chamber	kg	300			
	Max. loading per slide-in drip tray	kg	4	8		
	Max. loading per bottom drip tray	kg	4	8		
Textured stainless steel exterior	Width	(D) mm	824	1224		
	Height (sizes 450, 750, 1060 with castors)	(E) mm	1183	1247	1720	
	Depth (without door handle, depth of handle +56 mm)	(F) mm	684	784	1035	
Standard equipment	Standard works calibration certificate (measuring point chamber center)	°C	+37			
	CO <sub>2</sub> connection set: hose with coupling and clamp	number	2			
	Door-open-recognition		-	●		
Temperature	Working temperature range	°C	at least 10 above ambient temperature up to +80			
	Setting temperature range	°C	+20 to +80			
	Setting accuracy	°C	up to 99.9: 0.1 / from 100: 0.5			
Further data	Electrical load at 230 V, 50/60 Hz	W	840	-		
	Electrical load at 115 V, 50/60 Hz	W	840	-		
	Electrical load at 400 V and 3 x 230V w/o neutral, 50/60 Hz	W	-	1450	1725	
Packing data	Net weight	kg	110	161	217	260
	Gross weight (packed in carton)	kg	161	227	288	424
	Width	mm	930	1330		1370
	Height	mm	1380	1440	1910	1970
	Depth	mm	930	1050		1300

**Order No. Universal Ovens**

- U = Universal Oven
- N = Natural convection
- F = Forced air circulation
- plus = Model with TwinDISPLAY

UF260Pplus    UF450Pplus    UF750Pplus    UF1060Pplus




Options	260	450	750	1060
Voltage 115 V, 50/60 Hz	X2	-		
Full-sight glass door (4-layer insulating glass)				B0
Interior nearly gastight with possibility for gas inlet/outlet through 2 tubes with ball valves - technical clarification required for combination with other options				K3
Entry port with 23 mm clear diameter, can be closed, standard position	left centre/centre	F0		
	left centre/top	F1		
	right centre/centre	F2		
	right centre/top	F3		
Entry port with 23 mm clear diameter, can be closed, individual position (please state position)	left	F4		
	right	F5		
	rear	F6		
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please state location)				D6
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please state location)				F7
Entry port with 57 mm clear diameter, can be closed by flap, individual positions in the back wall (please state position, not possible for devices with forced air circulation of sizes 30-75)				F8
Entry port with 100 mm clear diameter, can be closed by flap, individual positions in the back wall (please state position, not possible for devices with forced air circulation of size 75)				F9
4 - 20 mA current loop interface (0 to +90 °C = 4 - 20 mA) Temperature controller, actual value Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 1 SingleDISPLAY, max. 3 TwinDISPLAY) - price per sensor				V3
				V6
Fan speed monitoring with switching off the heating and with alarm in case of failure (for UFplus)				V4
Door with lock and key (safety lock)				B6
Door hinged on the left	B8	-		
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)				H6
Potential-free contact (24 V/2 A) with socket, for signal generation, controlled by programme segment, for free-selectable functions to be activated (e.g. activation of audible and visual signals, exhaust motors, fans, stirrers, etc.). Only for units with TwinDISPLAY	2 contacts		H72	
Process-dependent programmable door lock (only for units with TwinDISPLAY)				D4
Flexible Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 1 additional sensor is possible). The measured temperature will be indicated on the display, recorded in the integral data store, and can be documented via the AtmoCONTROL software				H8
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)				C3
Castor frame (2-part), height 140 mm	R9	-		
Fresh-air filter (filtration efficiency 80 %) mounted at the device bottom (for devices with forced air circulation, sizes 30 - 260 require castor frame R9 or subframe)				R8

Accessories	260	450	750	1060
Reinforced stainless steel grid, electropolished, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber	B42444	B42447	B42450	
Reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber	B42445	B42448	B42451	
Stainless steel grid, electropolished	E28891	-	B41251	
Reinforced stainless steel grid, electropolished, max. loading 60 kg; from size 450 with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber	E29766	B32190	B32550	
Perforated stainless steel shelf	B29725	-		
Reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber	-	B32191	B32549	
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (may affect the temperature distribution, requires option K1)	-	B32763	-	
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1 for sizes 450/750)	B29726	-	B32599	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, requires option K1)	-	B34055	-	
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1 for sizes 450/750)	B29722	-	B29769	
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), straight, for exhaust air ducting (if necessary for connection by hose)				B29718

Accessories	260	450	750	1060
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose)	B29719			
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29738	B29740	B29742	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29739	B29741	B29743	
Set of height adjustable feet (4 pcs)	B29768	-		
Subframe, adjustable in height (size 30 to 75: height 600 mm, from size 110: height 500 mm)	B29751	B29753	-	
Subframe, adjustable in height, height 130 mm	B33664	-		

Accessories   Software and Connectivity	260	450	750	1060
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1			
Integration of one additional unit (up to max. 31 units) into an already existent FDA-software licence (only for units with TwinDISPLAY)	FDAQ2			
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)	B33170			
Ethernet connection cable 5 m for computer interface	E06189			
USB-Ethernet adapter	E06192			



**m360**

**MODIFICATIONS**

CONSTANT CLIMATE CHAMBERS HPPeco  
HUMIDITY CHAMBERS HCP  
CLIMATE CHAMBERS ICH  
ENVIRONMENTAL TEST CHAMBERS CTC / TTC  
UNIVERSAL OVENS U  
VACUUM OVENS VO  
INCUBATORS I  
CO<sub>2</sub> INCUBATORS ICO  
COMPRESSOR-COOLED INCUBATORS ICP  
PELTIER-COOLED INCUBATORS IPPeco



**CONSTANT CLIMATE CHAMBERS HPPeco**

Options   Mechanics	110	260	410	750	1060	1400	2200
Door handlebar   Reinforced ; price per door	-				Q4m		
Door hinge   Exchange		B8m			-		
Door hinge   Max. opening angle 120° ; price per door				Q1m			
Door hinge   Max. opening angle 95° ; price per door				Q0m			
Entry port   Inner diameter 14mm	back panel			N0m			
	side panel			N1m			
Entry port   Inner diameter 23mm	back panel		-			N2m	
Entry port   Inner diameter 38mm	back panel			N4m			
	side panel			N5m			
Entry port   Inner diameter 57mm	back panel			N6m			
	side panel			N7m			
Entry port   Inner diameter 80mm	back panel			N8m			
	side panel			N9m			
Entry port   Inner diameter 100mm	back panel			N10m			
	side panel			N11m			
Entry port   Inner diameter 120mm	back panel			N12m			
	side panel			N13m			
Entry port   Inner diameter 40mm ; Silicone	back panel		-			N20m	
Entry port   Inner diameter 80mm ; Silicone	back panel			N21m			
Entry port   Inner diameter 120mm ; Silicone	back panel			N22m			
ESD modification			M3m			-	
Heavy-duty   Interior modification option K1	-		K1		-		
Heavy-duty frame   Interior fittings ; requires base frame at extra cost	-	M11m			-		
		-		M22m		-	
		-			M31m	-	
Heavy-duty frame   Base frame ; requires option R9 for appliances size 260 at extra cost	-	M10m		-			
		-		M20m		-	
		-			M30m	-	
Heavy-duty frame   Grid ; price per grating ; requires base frame and interior fittings at extra cost	-	M12m		-			
		-		M23m		-	
		-			M32m	-	
Levelling castors   Max. height: 102mm ; requires option R9 for appliances up to and including size 260 at extra cost				M4m			
Levelling feet   Max. height: 140mm ; requires option R9 for appliances up to and including size 260 at extra cost				M5m			
Maintenance flap   Electronics box		M1m			-		
Metal outer door   Inspection window ; heated	-			B3m		-	
Metal outer door   Inspection window ; non-heated			B0m		-		
Seal   Housing ; ESD protection ; temperature resistance: 180 °C				Q3m		-	
Seal   Viton ; housing ; temperature resistance: 250 °C				Q2m		-	
Telescopic slides   Grid ; price per grid	M6m			-			
Telescopic slides   Perforated sheet ; price per perforated sheet	M7m			-			
Options   Electronics	110	260	410	750	1060	1400	2200
Control unit   External				E4m			
Control box   External				E5m			
Door   Electric lock ; safety switch		E7m			-		
Exhaust air flap   Rear panel				E6m			
Interface 4-20mA   Humidity				V7m			

Options   Electronics		110	260	410	750	1060	1400	2200
Lighting   Grid ; LED light strips ; cool white 6,500 K	first grid	J17m			-	J17m		
	additional grid	J17mx			-	J17mx		
Lighting   Grid ; LED light strips ; cool white 6,500 K + warm white 2,700 K	first grid	J18m			-	J18m		
	additional grid	J18mx			-	J18mx		
Lighting   Grid ; LED light strips ; warm white 2,700 K	first grid	J19m			-	J19m	-	
	additional grid	J19mx			-	J19mx		
Lighting   Grid reinforced ; LED light strips ; cool white 6,500 K ; requires option K1 at extra charge, except from size 1060	first grid	-	J27m			-		
	additional grid	-	J27mx			-		
Lighting   Grid reinforced ; LED light strips ; cool white 6,500 K + warm white 2,700 K ; requires option K1 at extra charge, except from size 1060	first grid	-	J28m			-		
	additional grid	-	J28mx			-		
Lighting   Grid reinforced ; LED light strips ; warm white 2,700 K ; requires option K1 at extra charge, except from size 1060	first grid	-	J29m			-		
	additional grid	-	J29mx			-		
Lighting   Interior K1 modification ; LED light strips ; cool white 6,500 K ; requires option K1 at extra charge, except from size 1060	first grid	-	J07m			-		
	additional grid	-	J07mx			-		
Lighting   Interior K1 modification ; LED light strips ; cool white 6,500 K + warm white 2,700 K ; requires option K1 at extra charge, except from size 1060	first grid	-	J08m			-		
	additional grid	-	J08mx			-		
Lighting   Interior K1 modification ; LED light strips ; warm white 2,700 K ; requires option K1 at extra charge, except from size 1060	first grid	-	J09m			-		
	additional grid	-	J09mx			-		
4 - 20 mA current loop interface	Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 3) - price per sensor				V6			
Parameter monitoring   CO <sub>2</sub> ; temperature resistance: +70 °C					U2m			
Parameter monitoring   O <sub>2</sub> ; temperature resistance: +200 °C					U0m			
Potential-free contact   Message TARGET-value reached					H5m			
Potential-free contact   Programmable ; only for appliances with TwinDISPLAY	2 contacts (A, B)				H7m2			
	4 contacts (A, B, C, D)	-				H7m4		
PT100 sensor   Freely positionable ; Connection socket on the rear of the appliance					H4m			
PT100 sensor   Freely positionable	1 sensor				H8m			
	2 sensors				H8m2			
	3 sensors				H8m3			
Signal lamp   Control panel ; red, yellow, green ; requires option H6 and H72 at extra cost		E0m			-			
Signal lamp   Top cover ; 3 lights (red, yellow, green) ; requires option H6 and H72 at extra cost		E1m			-			
Signal lamp   Top cover ; one multi-coloured light ; requires option H6 and H72 at extra cost		E2m			-			
Socket outlet   Interior ; separately switchable ; temperature resistance: +60 °C ; requires option A8 at extra cost					R4m			
Socket outlet   Interior ; separately switchable ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost					R44m			
Socket outlet   Interior ; switchable via appliance ; temperature resistance: +60 °C ; requires option A8 at extra cost					R3m			
Socket outlet   Interior ; switchable via appliance ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost					R33m			
Accessories		110	260	410	750	1060	1400	2200
Connection cable   5m ; 230V					B58607			
Floor grid   Heavy-duty		E36692	E36704		B31972	E40438	-	
Grid   Reinforced ; max. load 60kg					-	B59004		
Perforated sheet					-	B59007		



Accessories	110	260	410	750	1060	1400	2200
Perforated sheet   Reinforced ; max. load 100kg ; requires option K1	-	B58460		B58461	B58462		-
Perforated sheet   Reinforced ; max. load 80kg	-	B57086		B52092			-
Plug   Ø23mm ; silicone ; temperature resistance: 200°C				E05387			
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C				E07479			
Plug   Ø57mm ; silicone ; temperature resistance: 200°C				E05389			
Plug   Ø100mm ; silicone ; temperature resistance: 200°C				E05386			
Plug   Ø14mm ; silicone sponge ; temperature resistance: 180°C				E05790			
Plug   Ø23mm ; silicone sponge ; temperature resistance: 180°C				E05361			
Plug   Ø37-40mm ; silicone sponge ; temperature resistance: 180°C				E05362			
Plug   Ø57mm ; silicone sponge ; temperature resistance: 180°C				E07865			
Plug   Ø80mm ; silicone sponge ; temperature resistance: 180°C				E05625			
Plug   Ø100mm ; silicone sponge ; temperature resistance: 180°C				E07866			
Plug   Ø120mm ; silicone sponge ; temperature resistance: 180°C				E05920			
Stacking frame	B57998	B58000			-		
Telescopic slide   Grid ; price per grid ; requires option K1	-	B43102		B43093	B43100		-
Telescopic slide   Perforated sheet ; price per perforated sheet ; requires option K1	-	B51618		B44476	B56888		-



## HUMIDITY CHAMBERS HCP

Options   Electronics	50	105	150	240
Control unit   External	E4m			
Accessories	50	105	150	240
Connection cable   5m ; 230V	B58607			
Plug   Ø23mm ; silicone ; temperature resistance: 200°C	E05387			
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C	E07479			
Plug   Ø57mm ; silicone ; temperature resistance: 200°C	E05389			
Plug   Ø100mm ; silicone ; temperature resistance: 200°C	E05386			
Plug   Ø14mm ; silicone sponge ; temperature resistance: 180°C	E05790			
Plug   Ø23mm ; silicone sponge ; temperature resistance: 180°C	E05361			
Plug   Ø37-40mm ; silicone sponge ; temperature resistance: 180°C	E05362			
Plug   Ø57mm ; silicone sponge ; temperature resistance: 180°C	E07865			
Plug   Ø80mm ; silicone sponge ; temperature resistance: 180°C	E05625			
Plug   Ø100mm ; silicone sponge ; temperature resistance: 180°C	E07866			
Plug   Ø120mm ; silicone sponge ; temperature resistance: 180°C	E05920			



## CLIMATE CHAMBERS ICH

Options   Mechanics		110	260	750
Door handlebar   Reinforced ; price per door		-		Q4m
Door hinge   Exchange		B8m		-
Door hinge   Max. opening angle 120° ; price per door			Q1m	
Door hinge   Max. opening angle 95° ; price per door			Q0m	
Entry port   Inner diameter 14mm	back panel		N0m	
	side panel		N1m	
Entry port   Inner diameter 23mm	back panel		N2m	
	side panel		N3m	
Entry port   Inner diameter 38mm	back panel		N4m	
	side panel		N5m	
Entry port   Inner diameter 57mm	back panel		N6m	
	side panel		N7m	
Entry port   Inner diameter 80mm	back panel		N8m	
	side panel		N9m	
Entry port   Inner diameter 100mm	back panel		N10m	
	side panel		N11m	
Entry port   Inner diameter 120mm	back panel		N12m	
	side panel		N13m	
Heavy-duty   Interior modification option K1		-	K1	-
Levelling castors   Max. height: 102mm ; requires option R9 for appliances up to and including size 260 at extra cost			M4m	
Levelling feet   Max. height: 140mm ; requires option R9 for appliances up to and including size 260 at extra cost			M5m	
Metal outer door   Inspection window ; heated		-		B3m
Metal outer door   Inspection window ; non-heated			B0m	
Seal   Housing ; ESD protection ; temperature resistance: 180 °C			Q3m	
Seal   Viton ; housing ; temperature resistance: 250 °C			Q2m	
Options   Electronics		110	260	750
Control unit   External			E4m	
Control box   External			E5m	
Door   Electric lock ; safety switch		E7m		-
Exhaust air flap   Rear panel			E6m	
Interface 4-20mA   Humidity			V7m	
Lighting   Interior			R0m	
Parameter monitoring   CO <sub>2</sub> ; temperature resistance: +70 °C			U2m	
Parameter monitoring   O <sub>2</sub> ; temperature resistance: +200 °C			U0m	
Potential-free contact   Message TARGET-value reached			H5m	
Potential-free contact   Programmable ; only for appliances with TwinDISPLAY	2 contacts (A, B)		H7m2	
	4 contacts (A, B, C, D)		-	H7m4
PT100 sensor   Freely positionable ; Connection socket on the rear of the appliance			H4m	
PT100 sensor   Freely positionable	1 sensor		H8m	
	2 sensors		H8m2	
	3 sensors		H8m3	
Signal lamp   Control panel ; red, yellow, green ; requires option H6 and H72 at extra cost			E0m	
Signal lamp   Top cover ; 3 lights (red, yellow, green) ; requires option H6 and H72 at extra cost			E1m	
Signal lamp   Top cover ; one multi-coloured light ; requires option H6 and H72 at extra cost			E2m	
Socket outlet   Interior ; separately switchable ; temperature resistance: +60 °C ; requires option A8 at extra cost			R4m	
Socket outlet   Interior ; separately switchable ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost			R44m	
Socket outlet   Interior ; switchable via appliance ; temperature resistance: +60 °C ; requires option A8 at extra cost			R3m	
Socket outlet   Interior ; switchable via appliance ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost			R33m	

Accessories	110	260	750
Connection cable   5m ; 230V		B58607	
Floor grid   Heavy-duty	E36692	E36704	B31972
Perforated sheet   Reinforced ; max. load 100kg ; requires option K1	-	B58460	B58461
Perforated sheet   Reinforced ; max. load 80kg	-	B57086	B52092
Plug   Ø23mm ; silicone ; temperature resistance: 200°C		E05387	
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C		E07479	
Plug   Ø57mm ; silicone ; temperature resistance: 200°C		E05389	
Plug   Ø100mm ; silicone ; temperature resistance: 200°C		E05386	
Plug   Ø14mm ; silicone sponge ; temperature resistance: 180°C		E05790	
Plug   Ø23mm ; silicone sponge ; temperature resistance: 180°C		E05361	
Plug   Ø37-40mm ; silicone sponge ; temperature resistance: 180°C		E05362	
Plug   Ø57mm ; silicone sponge ; temperature resistance: 180°C		E07865	
Plug   Ø80mm ; silicone sponge ; temperature resistance: 180°C		E05625	
Plug   Ø100mm ; silicone sponge ; temperature resistance: 180°C		E07866	
Plug   Ø120mm ; silicone sponge ; temperature resistance: 180°C		E05920	



## ENVIRONMENTAL TEST CHAMBERS CTC / TTC

Options   Electronics		CTC256	TTC256
PT100 sensor   Freely positionable ; Connection socket on the rear of the appliance		H4m	-
PT100 sensor   Freely positionable	1 sensor	H8m	-
Accessories		CTC256	TTC256
Floor grid   Heavy-duty		E41309	
Plug   Ø23mm ; silicone ; temperature resistance: 200°C		E05387	
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C		E07479	
Plug   Ø57mm ; silicone ; temperature resistance: 200°C		E05389	
Plug   Ø100mm ; silicone ; temperature resistance: 200°C		E05386	
Plug   Ø14mm ; silicone sponge ; temperature resistance: 180°C		E05790	
Plug   Ø23mm ; silicone sponge ; temperature resistance: 180°C		E05361	
Plug   Ø37-40mm ; silicone sponge ; temperature resistance: 180°C		E05362	
Plug   Ø57mm ; silicone sponge ; temperature resistance: 180°C		E07865	
Plug   Ø80mm ; silicone sponge ; temperature resistance: 180°C		E05625	
Plug   Ø100mm ; silicone sponge ; temperature resistance: 180°C		E07866	
Plug   Ø120mm ; silicone sponge ; temperature resistance: 180°C		E05920	

## UNIVERSAL OVENS U

Options   Mechanics	30	55	75	110	160	260	450	750	1060
Door handlebar   Reinforced ; price per door	-								Q4m
Door hinge   Max. opening angle 120° ; price per door	-					Q1m			
Door hinge   Max. opening angle 95° ; price per door	-					Q0m			
Door lock   Magnetic ; temperature resistance: 200 °C ; requires option A8 at extra cost	-					M0m		-	
Entry port   Inner diameter 14mm	side panel		N51m						
Entry port   Inner diameter 38mm	side panel		N56m						
Entry port   Inner diameter 57mm	back panel		F8m		-				
	side panel		N59m						
Entry port   Inner diameter 80mm	back panel		N60m						
	side panel		N61m						
Entry port   Inner diameter 100mm	back panel		F9m		-				
	side panel		N63m						
Entry port   Inner diameter 120mm	back panel		N64m						
	side panel		N65m						
ESD modification	-		M3m						
Filter holder (one filter)   Air inlet ; filter available separately as an accessory ; requires a base frame or option R9 for appliances up to and including size 260	-					R81m			
Filter holder (two filters)   Air inlet ; filter available separately as an accessory ; requires a base frame or option R9 for appliances up to and including size 260	-					R82m			
Heavy-duty   Interior modification option K1	-					K1		-	
Heavy-duty frame   Interior fittings ; requires base frame at extra cost	-					M11m		-	
	-					M21m			-
	-					M22m			-
	-					M31m			
Heavy-duty frame   Base frame ; requires option R9 for appliances size 260 at extra cost	-					M10m		-	
	-					M20m			-
	-					M30m			
Heavy-duty frame   Grid ; price per grating ; requires base frame and interior fittings at extra cost	-					M12m		-	
	-					M23m			-
	-					M32m			
Insulation   Modification for particle reduction	-					M2m			
Levelling castors   Max. height: 102mm ; requires option R9 for appliances up to and including size 260 at extra cost	-					M4m			
Levelling feet   Max. height: 140mm ; requires option R9 for appliances up to and including size 260 at extra cost	-					M5m			
Maintenance flap   Electronics box	-		M1m			-			
Seal   Housing ; ESD protection ; temperature resistance: 180 °C ; requires option A8 at extra cost	-					Q3m			
Seal   Viton ; housing ; temperature resistance: 250 °C ; requires option A8 at extra cost	-					Q2m			
Telescopic slides   Grid ; price per grid	-		M6m			-			
Telescopic slides   Perforated sheet ; price per perforated sheet	-		M7m			-			
Options   Electronics	30	55	75	110	160	260	450	750	1060
Control unit   External	-					E4m			
Control box   External	-		E5m						
Door   Electric lock ; safety switch	-		E7m			-			



Options   Electronics	30	55	75	110	160	260	450	750	1060	
Interface 4-20mA   Humidity ; requires option U1m at extra cost	V7m									
Exhaust air extraction   Rapid cooling fan	-		E3m							
Parameter monitoring   CO <sub>2</sub> ; temperature resistance: +70 °C ; requires option A8 at extra cost	U2m									
Parameter monitoring   Humidity ; temperature resistance +180°C ; requires option A8 at extra cost	U1m									
Parameter monitoring   O <sub>2</sub> ; temperature resistance: +200 °C ; requires option A8 at extra cost	U0m									
Potential-free contact   Message TARGET-value reached	H5m									
Potential-free contact   Programmable ; only for appliances with Twin-DISPLAY	4 contacts (A, B, C, D)						-			
PT100 sensor   Freely positionable ; Connection socket on the rear of the appliance	H4m									
PT100 sensor   Freely positionable	2 sensors		H8m2							
	3 sensors		H8m3							
Signal lamp   Control panel ; red, yellow, green ; requires option H6 and H72 at extra cost	-			E0m						
Signal lamp   Top cover ; 3 lights (red, yellow, green) ; requires option H6 and H72 at extra cost	-			E1m						
Signal lamp   Top cover ; one multi-coloured light ; requires option H6 and H72 at extra cost	-			E2m						
Socket outlet   Interior ; separately switchable ; temperature resistance: +60 °C ; requires option A8 at extra cost	R4m									
Socket outlet   Interior ; separately switchable ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost	R44m									
Socket outlet   Interior ; switchable via appliance ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost	R33m									
Accessories	30	55	75	110	160	260	450	750	1060	
Air filter   Supply air opening ; filter class: F7 ; requires option R81m / R82m	E07438									
Air filter   Supply air opening ; filter class: H13 ; requires option R81m / R82m	E07439									
Connection cable   5m ; 230V	B58607						-			
Exhaust air   Rear wall ; condensate drainage	B37684									
Floor grid   Heavy-duty	E50944	E46407	E36692		E36704	B31972	E40438			
Perforated sheet   Reinforced ; max. load 100kg ; requires option K1	-					B58460	B58461	B58462		
Perforated sheet   Reinforced ; max. load 80kg	-					B57086	B52092	-		
Plug   Ø23mm ; silicone ; temperature resistance: 200°C	E05387									
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C	E07479									
Plug   Ø37-40mm ; silicone ; temperature resistance: 300°C	E07699									
Plug   Ø57mm ; silicone ; temperature resistance: 200°C	E05389									
Plug   Ø100mm ; silicone ; temperature resistance: 200°C	E05386									
Plug   Ø23mm ; silicone ; temperature resistance: 300°C	E07698									
Plug   Ø57mm ; silicone ; temperature resistance: 300°C	E07700									
Plug   Ø80mm ; silicone ; temperature resistance: 300°C	E07701									
Plug   Ø100-120mm ; silicone ; temperature resistance: 300°C	E07702									

Accessories	30	55	75	110	160	260	450	750	1060
Stacking frame	-		B50598	B50599		B50600	-		
Telescopic slide   Grid ; price per grid ; requires option K1			-			B43102	B43093	B43100	
Telescopic slide   Perforated sheet ; price per perforated sheet ; requires option K1			-			B51618	B44476	B56888	



VACUUM OVENS VO				
Options   Mechanics		29	49	101
Door hinge   Exchange		B8m		
Door lock   Magnetic ; temperature resistance: 200 °C		M0m		
Entry port   DN16 ISO-KF		N30m		
Entry port   DN25 ISO-KF		N31m		
Entry port   DN40 ISO-KF		N32m		
Entry port   DN50 ISO-KF		N33m		
Options   Electronics		29	49	101
Control unit   External		E4m		
Door   Opening detection		V5m		
Lighting   Interior		R0m		
Potential-free contact   Programmable ; only for appliances with TwinDISPLAY		1 contact (A)		H7m
		2 contacts (A, B)		H7m2
		3 contacts (A, B, C)		H7m3
PT100 sensor   Freely positionable		H8m		
Signal lamp   Control panel ; red, yellow, green ; requires option H6 and H72 at extra cost		-	E0m	
Signal lamp   Top cover ; 3 lights (red, yellow, green) ; requires option H6 and H72 at extra cost		-	E1m	
Signal lamp   Top cover ; one multi-coloured light ; requires option H6 and H72 at extra cost		E2m		
Accessories		29	49	101
Connection cable   5m ; 230V		B58607		
Stapelgestell   Nivellierbar		B50577	B59029	



## INCUBATORS I

Options   Mechanics	30	55	75	110	160	260	450	750
Door handlebar   Reinforced ; price per door								Q4m
Door hinge   Max. opening angle 120° ; price per door				Q1m				
Door hinge   Max. opening angle 95° ; price per door				Q0m				
Door lock   Magnetic ; temperature resistance: 200 °C			M0m					-
Entry port   Inner diameter 14mm side panel				N51m				
Entry port   Inner diameter 38mm side panel				N56m				
Entry port   Inner diameter 57mm back panel		F8m				-		
				N59m				
Entry port   Inner diameter 80mm back panel				N60m				
				N61m				
Entry port   Inner diameter 100mm back panel		F9m				-		
				N63m				
Entry port   Inner diameter 120mm back panel				N64m				
				N65m				
ESD modification		-				M3m		
Filter holder (one filter)   Air inlet ; filter available separately as an accessory ; requires a base frame or option R9 for appliances up to and including size 260				R81m				
Filter holder (two filters)   Air inlet ; filter available separately as an accessory ; requires a base frame or option R9 for appliances up to and including size 260				R82m				
Heavy-duty   Interior modification option K1			-			K1		-
Heavy-duty frame   Interior fittings ; requires base frame at extra cost			-			M11m		-
			-				M21m	-
			-					M22m
Heavy-duty frame   Base frame ; requires option R9 for appliances size 260 at extra cost	-	M20m		-				M20m
						M10m		-
Heavy-duty frame   Grid ; price per grating ; requires base frame and interior fittings at extra cost			-			M12m		-
			-					M23m
Insulation   Modification for particle reduction				M2m				
Levelling castors   Max. height: 102mm ; requires option R9 for appliances up to and including size 260 at extra cost				M4m				
Levelling feet   Max. height: 140mm ; requires option R9 for appliances up to and including size 260 at extra cost				M5m				
Maintenance flap   Electronics box		-			M1m			-
Seal   Housing ; ESD protection ; temperature resistance: 180 °C				Q3m				
Seal   Viton ; housing ; temperature resistance: 250 °C				Q2m				
Telescopic slides   Grid ; price per grid	-		M6m					-
Telescopic slides   Perforated sheet ; price per perforated sheet	-		M7m					-
Options   Electronics	30	55	75	110	160	260	450	750
Control unit   External				E4m				
Control box   External		-				E5m		
Door   Electric lock ; safety switch		-			E7m			-
Interface 4-20mA   Humidity ; requires option U1m at extra cost				V7m				
Exhaust air extraction   Rapid cooling fan		-				E3m		
Parameter monitoring   CO <sub>2</sub> ; temperature resistance: +70 °C				U2m				
Parameter monitoring   Humidity ; temperature resistance +180°C				U1m				
Parameter monitoring   O <sub>2</sub> ; temperature resistance: +200 °C				U0m				
Potential-free contact   Message TARGET-value reached				H5m				



Options   Electronics	30	55	75	110	160	260	450	750
Potential-free contact   Programmable ; only for appliances with TwinDISPLAY	2 contacts (A, B)		H7m2					
	4 contacts (A, B, C, D)		-				H7m4	
PT100 sensor   Freely positionable ; Connection socket on the rear of the appliance	H4m							
PT100 sensor   Freely positionable	2 sensors		H8m2					
	3 sensors		H8m3					
Signal lamp   Control panel ; red, yellow, green ; requires option H6 and H72 at extra cost	-		E0m					
Signal lamp   Top cover ; 3 lights (red, yellow, green) ; requires option H6 and H72 at extra cost	-		E1m					
Signal lamp   Top cover ; one multi-coloured light ; requires option H6 and H72 at extra cost	-		E2m					
Socket outlet   Interior ; separately switchable ; temperature resistance: +60 °C ; requires option A8 at extra cost	R4m							
Socket outlet   Interior ; separately switchable ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost	R44m							
Socket outlet   Interior ; switchable via appliance ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost	R33m							
Accessories	30	55	75	110	160	260	450	750
Air filter   Supply air opening ; filter class: F7 ; requires option R81m / R82m	E07438							
Air filter   Supply air opening ; filter class: H13 ; requires option R81m / R82m	E07439							
Connection cable   5m ; 230V	B58607							
Exhaust air   Rear wall ; condensate drainage	B37684							
Floor grid   Heavy-duty	E50944	E46407		E36692		E36704	B31972	
Perforated sheet   Reinforced ; max. load 100kg ; requires option K1	-		-			B58460	B58461	
Perforated sheet   Reinforced ; max. load 80kg	-		-			B57086	B52092	
Plug   Ø23mm ; silicone ; temperature resistance: 200°C	E05387							
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C	E07479							
Plug   Ø57mm ; silicone ; temperature resistance: 200°C	E05389							
Plug   Ø100mm ; silicone ; temperature resistance: 200°C	E05386							
Plug   Ø14mm ; silicone sponge ; temperature resistance: 180°C	E05790							
Plug   Ø23mm ; silicone sponge ; temperature resistance: 180°C	E05361							
Plug   Ø37-40mm ; silicone sponge ; temperature resistance: 180°C	E05362							
Plug   Ø57mm ; silicone sponge ; temperature resistance: 180°C	E07865							
Plug   Ø80mm ; silicone sponge ; temperature resistance: 180°C	E05625							
Plug   Ø100mm ; silicone sponge ; temperature resistance: 180°C	E07866							
Plug   Ø120mm ; silicone sponge ; temperature resistance: 180°C	E05920							
Stacking frame	-		B50598	B50599		B50600	-	
Telescopic slide   Grid ; price per grid ; requires option K1	-		-			B43102	B43093	
Telescopic slide   Perforated sheet ; price per perforated sheet ; requires option K1	-		-			B51618	B44476	

CO<sub>2</sub> INCUBATORS ICO

Options   Electronics	50	105	150	240
Control unit   External			E4m	
Door   Electric lock			D4m	
Accessories	50	105	150	240
Connection cable   5m ; 230V			B58607	
Plug   Ø23mm ; silicone ; temperature resistance: 200°C			E05387	
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C			E07479	
Plug   Ø57mm ; silicone ; temperature resistance: 200°C			E05389	
Plug   Ø100mm ; silicone ; temperature resistance: 200°C			E05386	
Plug   Ø14mm ; silicone sponge ; temperature resistance: 180°C			E05790	
Plug   Ø23mm ; silicone sponge ; temperature resistance: 180°C			E05361	
Plug   Ø37-40mm ; silicone sponge ; temperature resistance: 180°C			E05362	
Plug   Ø57mm ; silicone sponge ; temperature resistance: 180°C			E07865	
Plug   Ø80mm ; silicone sponge ; temperature resistance: 180°C			E05625	
Plug   Ø100mm ; silicone sponge ; temperature resistance: 180°C			E07866	
Plug   Ø120mm ; silicone sponge ; temperature resistance: 180°C			E05920	



## COMPRESSOR-COOLED INCUBATORS ICP

Options   Mechanics		110	260	450	750
Door handlebar   Reinforced ; price per door			-		Q4m
Door hinge   Exchange		B8m		-	
Door hinge   Max. opening angle 120° ; price per door			Q1m		
Door hinge   Max. opening angle 95° ; price per door			Q0m		
Entry port   Inner diameter 14mm	back panel		N0m		
	side panel		N1m		
Entry port   Inner diameter 23mm	back panel		N2m		
	side panel		N3m		
Entry port   Inner diameter 38mm	back panel		N4m		
	side panel		N5m		
Entry port   Inner diameter 57mm	back panel		N6m		
	side panel		N7m		
Entry port   Inner diameter 80mm	back panel		N8m		
	side panel		N9m		
Entry port   Inner diameter 100mm	back panel		N10m		
	side panel		N11m		
Entry port   Inner diameter 120mm	back panel		N12m		
	side panel		N13m		
Heavy-duty   Interior modification option K1		-	K1		-
Levelling castors   Max. height: 102mm ; requires option R9 for appliances up to and including size 260 at extra cost			M4m		
Levelling feet   Max. height: 140mm ; requires option R9 for appliances up to and including size 260 at extra cost			M5m		
Metal outer door   Inspection window ; heated			-		B3m
Metal outer door   Inspection window ; non-heated		B0m		-	B0m
Seal   Housing ; ESD protection ; temperature resistance: 180 °C			Q3m		
Seal   Viton ; housing ; temperature resistance: 250 °C			Q2m		
Options   Electronics		110	260	450	750
Control unit   External			E4m		
Control box   External			E5m		
Door   Electric lock ; safety switch		E7m		-	
Exhaust air flap   Rear panel			E6m		
Interface 4-20mA   Humidity ; requires option U1m at extra cost			V7m		
Compressed air drying   Dehumidification module ; requires option U1m at extra cost			C9m		
Lighting   Interior			R0m		
Parameter monitoring   CO <sub>2</sub> ; temperature resistance: +70 °C			U2m		
Parameter monitoring   Humidity ; temperature resistance +180°C			U1m		
Potential-free contact   Message TARGET-value reached			H5m		
Potential-free contact   Programmable ; only for appliances with TwinDISPLAY	2 contacts (A, B)		H7m2		
	4 contacts (A, B, C, D)		-	H7m4	
PT100 sensor   Freely positionable ; Connection socket on the rear of the appliance			H4m		
PT100 sensor   Freely positionable	1 sensor		H8m		
	2 sensors		H8m2		
	3 sensors		H8m3		
Signal lamp   Control panel ; red, yellow, green ; requires option H6 and H72 at extra cost			E0m		
Signal lamp   Top cover ; 3 lights (red, yellow, green) ; requires option H6 and H72 at extra cost			E1m		
Signal lamp   Top cover ; one multi-coloured light ; requires option H6 and H72 at extra cost			E2m		
Socket outlet   Interior ; separately switchable ; temperature resistance: +60 °C ; requires option A8 at extra cost			R4m		
Socket outlet   Interior ; separately switchable ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost			R44m		

Options   Electronics	110	260	450	750
Socket outlet   Interior ; switchable via appliance ; temperature resistance: +60 °C ; requires option A8 at extra cost	R3m			
Socket outlet   Interior ; switchable via appliance ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost	R33m			
Accessories	110	260	450	750
Connection cable   5m ; 230V	B58607			
Floor grid   Heavy-duty	E36692	E36704	B31972	
Perforated sheet   Reinforced ; max. load 100kg ; requires option K1	-	B58460	B58461	
Perforated sheet   Reinforced ; max. load 80kg	-	B57086	B52092	
Plug   Ø23mm ; silicone ; temperature resistance: 200°C	E05387			
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C	E07479			
Plug   Ø57mm ; silicone ; temperature resistance: 200°C	E05389			
Plug   Ø100mm ; silicone ; temperature resistance: 200°C	E05386			
Plug   Ø14mm ; silicone sponge ; temperature resistance: 180°C	E05790			
Plug   Ø23mm ; silicone sponge ; temperature resistance: 180°C	E05361			
Plug   Ø37-40mm ; silicone sponge ; temperature resistance: 180°C	E05362			
Plug   Ø57mm ; silicone sponge ; temperature resistance: 180°C	E07865			
Plug   Ø80mm ; silicone sponge ; temperature resistance: 180°C	E05625			
Plug   Ø100mm ; silicone sponge ; temperature resistance: 180°C	E07866			
Plug   Ø120mm ; silicone sponge ; temperature resistance: 180°C	E05920			



**PELTIER-COOLED INCUBATORS IPPeco**

Options   Mechanics	30	55	110	260	410	750	1060	1400	2200
Door handlebar   Reinforced ; price per door	-				Q4m				
Door hinge   Exchange	B8m				-				
Door hinge   Max. opening angle 120° ; price per door					Q1m				
Door hinge   Max. opening angle 95° ; price per door					Q0m				
Entry port   Inner diameter 14mm	back panel					D6m			
	side panel					N51m			
Entry port   Inner diameter 23mm	back panel	-				F6m			
	side panel					N56m			
Entry port   Inner diameter 38mm	back panel	-				F7m			
	side panel					N59m			
Entry port   Inner diameter 57mm	back panel					F8m			
	side panel					N60m			
Entry port   Inner diameter 80mm	back panel					N61m			
	side panel					F9m			
Entry port   Inner diameter 100mm	back panel					N63m			
	side panel					F9m			
Entry port   Inner diameter 40mm ; Silicone	back panel	-				N20m			
Entry port   Inner diameter 80mm ; Silicone	back panel	-				N21m			
Entry port   Inner diameter 120mm ; Silicone	back panel					N22m			
Entry port   Inner diameter 120mm	back panel					N64m			
	side panel					N65m			
ESD modification	-					M3m			-
Heavy-duty   Interior modification option K1	-			K1					-
Heavy-duty frame   Interior fittings ; requires base frame at extra cost	-		M11m					-	
					M22m	-			
							M31m	-	
Heavy-duty frame   Base frame ; requires option R9 for appliances size 260 at extra cost	-		M10m					-	
					M20m	-			
							M30m	-	
Heavy-duty frame   Grid ; price per grating ; requires base frame and interior fittings at extra cost	-		M12m					-	
					M23m	-			
							M32m	-	
Levelling castors   Max. height: 102mm ; requires option R9 for appliances up to and including size 260 at extra cost					M4m				
Levelling feet   Max. height: 140mm ; requires option R9 for appliances up to and including size 260 at extra cost					M5m				
Maintenance flap   Electronics box	-	M1m							-
Metal outer door   Inspection window ; heated	-				B3m			-	
Metal outer door   Inspection window ; non-heated					B0m				
Seal   Housing ; ESD protection ; temperature resistance: 180 °C					Q3m				
Seal   Viton ; housing ; temperature resistance: 250 °C					Q2m				
Telescopic slides   Grid ; price per grid	-	M6m						-	
Telescopic slides   Perforated sheet ; price per perforated sheet	-	M7m						-	
Options   Electronics	30	55	110	260	410	750	1060	1400	2200
Control unit   External					E4m				
Control box   External	-					E5m			

Options   Electronics	30	55	110	260	410	750	1060	1400	2200
Door   Electric lock ; safety switch	-					E7m			
Exhaust air flap   Rear panel						E6m			
Interface 4-20mA   Humidity ; requires option U1m at extra cost	-					V7m			
Lighting   Grid ; LED light strips ; cool white 6,500 K	first grid	-			J17m		-	J17m	
	additional grid	-			J17mx		-	J17mx	
Compressed air drying   Dehumidification module ; requires option U1m at extra cost				C9m					-
Lighting   Grid ; LED light strips ; cool white 6,500 K + warm white 2,700 K	first grid	-			J18m		-	J18m	
	additional grid	-			J18mx		-	J18mx	
Lighting   Grid ; LED light strips ; warm white 2,700 K	first grid	-			J19m		-	J19m	-
	additional grid	-			J19mx		-	J19mx	
Lighting   Grid reinforced ; LED light strips ; cool white 6,500 K ; requires option K1 at extra charge, except from size 1060	first grid		-			J27m			-
	additional grid		-			J27mx			-
Lighting   Grid reinforced ; LED light strips ; cool white 6,500 K + warm white 2,700 K ; requires option K1 at extra charge, except from size 1060	first grid		-			J28m			-
	additional grid		-			J28mx			-
Lighting   Grid reinforced ; LED light strips ; warm white 2,700 K ; requires option K1 at extra charge, except from size 1060	first grid		-			J29m			-
	additional grid		-			J29mx			-
Lighting   Interior K1 modification ; LED light strips ; cool white 6,500 K ; requires option K1 at extra charge, except from size 1060	first grid		-			J07m			
	additional grid		-			J07mx			
Lighting   Interior K1 modification ; LED light strips ; cool white 6,500 K + warm white 2,700 K ; requires option K1 at extra charge, except from size 1060	first grid		-			J08m			
	additional grid		-			J08mx			
Lighting   Interior K1 modification ; LED light strips ; warm white 2,700 K ; requires option K1 at extra charge, except from size 1060	first grid		-			J09m			
	additional grid		-			J09mx			
4 - 20 mA current loop interface						V6			
Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 1 SingleDISPLAY, max. 3 TwinDISPLAY) - price per sensor									
Heavy-duty   Interior modification option K1		-			K1				-
Parameter monitoring   CO <sub>2</sub> ; temperature resistance: +70 °C						U2m			
Parameter monitoring   Humidity ; temperature resistance +180°C						U1m			
Parameter monitoring   O <sub>2</sub> ; temperature resistance: +200 °C						U0m			
Potential-free contact   Message TARGET-value reached						H5m			
Potential-free contact   Programmable ; only for appliances with Twin-DISPLAY	2 contacts (A, B)					H7m2			
	1 contact (A)					H7m			
	4 contacts (A, B, C, D)		-			H7m4			
PT100 sensor   Freely positionable ; Connection socket on the rear of the appliance						H4m			



Options   Electronics	30	55	110	260	410	750	1060	1400	2200
PT100 sensor   Freely positionable	1 sensor				H8m				
	2 sensors				H8m2				
	3 sensors				H8m3				
Signal lamp   Control panel ; red, yellow, green ; requires option H6 and H72 at extra cost					E0m				
Signal lamp   Top cover ; 3 lights (red, yellow, green) ; requires option H6 and H72 at extra cost					E1m				
Signal lamp   Top cover ; one multi-coloured light ; requires option H6 and H72 at extra cost					E2m				
Socket outlet   Interior ; separately switchable ; temperature resistance: +60 °C ; requires option A8 at extra cost					R4m				
Socket outlet   Interior ; separately switchable ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost					R44m				
Socket outlet   Interior ; switchable via appliance ; temperature resistance: +60 °C ; requires option A8 at extra cost					R3m				
Socket outlet   Interior ; switchable via appliance ; UK ; temperature resistance: +40 °C ; requires option A8 at extra cost					R33m				

Accessories	30	55	110	260	410	750	1060	1400	2200	
Connection cable   5m ; 230V	B58607									
Floor grid   Heavy-duty	E50944	E46407	E36692	E36704		B31972	E40438	-		
Grid   Reinforced ; max. load 60kg									B59004	
Perforated sheet									B59007	
Perforated sheet   Reinforced ; max. load 100kg ; requires option K1	-				B58460	B58461	B58462	-		
Perforated sheet   Reinforced ; max. load 80kg	-				B57086	B52092	-			
Plug   Ø23mm ; silicone ; temperature resistance: 200°C					E05387					
Plug   Ø37-40mm ; silicone ; temperature resistance: 200°C					E07479					
Plug   Ø57mm ; silicone ; temperature resistance: 200°C					E05389					
Plug   Ø100mm ; silicone ; temperature resistance: 200°C					E05386					
Plug   Ø14mm ; silicone sponge ; temperature resistance: 180°C					E05790					
Plug   Ø23mm ; silicone sponge ; temperature resistance: 180°C					E05361					
Plug   Ø37-40mm ; silicone sponge ; temperature resistance: 180°C					E05362					
Plug   Ø57mm ; silicone sponge ; temperature resistance: 180°C					E07865					
Plug   Ø80mm ; silicone sponge ; temperature resistance: 180°C					E05625					
Plug   Ø100mm ; silicone sponge ; temperature resistance: 180°C					E07866					
Plug   Ø120mm ; silicone sponge ; temperature resistance: 180°C					E05920					
Stacking frame	-		B57998	B58000					-	
Telescopic slide   Grid ; price per grid ; requires option K1	-			B43102		B43093	B43100	-		
Telescopic slide   Perforated sheet ; price per perforated sheet ; requires option K1	-			B51618		B44476	B56888	-		



PERSONAL NOTES

A large rectangular area filled with horizontal dotted lines, intended for taking personal notes. The lines are evenly spaced and extend across the width of the page.




**TECHNICAL  
INFORMATIONS  
MECHANICS**

Setup | Levelling castors  
Setup | Levelling feet  
Entry port | Moisture tight  
Entry port | KF  
Entry port | Silicone  
Entry ports  
Modification | ESD protection  
Modification | Insulation protection  
Heavy-duty | Base frame  
Heavy-duty | Grid (for rack)  
Heavy-duty | Interior frame  
Telescopic slides | Sizes 55 to 160  
Door | Limit opening angle  
Door | Replace door stop  
Door | Seals  
Door | Handlebar  
Door | Magnetic door lock  
Door | Full-sight glass  
Maintenance | Electronics box  
Supply air | Filter system



## Setup | Levelling castors


On request, the castors on Memmert appliances can be replaced with heavy-duty levelling castors. This allows the appliance to be moved to the desired installation location. On uneven surfaces, the levelling castors can be adjusted to level the floor.



Model sizes/Description		Levelling castors
Description	Specification	Castor with levelling feet
Technical data	Technical feature	movable
	Maximum height	mm 102
	Minimum height	mm 91
Ordering information	Further options (only after technical clarification)	ESD protection design or cleanroom suitability after technical clarification
	Retrofitting (after technical clarification)	Retrofitting possible
	Required option (depending on device type)	R9
<b>Order no.</b>		M4m

### Setup | Levelling feet

On request, the castors on Memmert appliances can be replaced with levelling feet. This ensures that users cannot move the appliance. On uneven surfaces, the levelling feet can be adjusted to compensate the unevenness in the floor.




Model sizes/Description		Levelling feet
Description	Specification	Levelling feet
Technical data	Technical feature	not movable
	Maximum height	mm
	Minimum height	mm
Ordering information	Further options (only after technical clarification)	<ul style="list-style-type: none"> <li>Cleanroom suitability possible after technical clarification.</li> <li>Customised heights on request</li> </ul>
	Retrofitting (after technical clarification)	Retrofitting possible
	Required option (depending on device type)	R9
<b>Order no.</b>		M5m



## Entry port | Moisture tight


Moisture-tight entry ports are mainly used, for example, to insert sensors into the appliance. This entry port is specially welded to prevent moisture escaping from the chamber and to ensure that the climate remains constant.



Model sizes/Description			Ø14	Ø23	Ø38	Ø57	Ø80	Ø100	Ø120						
Description	Specification	Entry ports for appliances with humidity control or special sealing requirements													
Technical data	Clear diameter	mm	14	23	38	57	80	100	120						
	Technical feature	moisture tight													
	Side panel	-	●	-	●	-	●	-	●	-	●				
	Back panel	●	-	●	-	●	-	●	-	●	-				
	Closure	Flap incl. plug					Slider incl. plug								
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible													
<b>Order no.</b>		N0m	N1m	N2m	N3m	N4m	N5m	N6m	N7m	N8m	N9m	N10m	N11m	N12m	N13m

### Entry port | KF

These entry ports were specially designed for vacuum ovens. The entry ports are welded vacuum-tight. This means that measuring sensors can be inserted into the device without pressure loss during vacuum operation.




Model sizes/Description		DN16	DN25	DN40	DN50
Description	Specification	DN16	DN25	DN40	DN50
	Design	ISO-KF			
Technical data	Technical feature	pressure tight			
	Side panel	-			
	Back panel	●			
	Closure	ISO-KF small flange			
Ordering information	Further options (only after technical clarification)	<ul style="list-style-type: none"> <li>• ISO-CF connection</li> <li>• further pressure-resistant connections, e.g. power supply</li> <li>• further sizes on request</li> <li>• further product series on request</li> <li>• Lateral insertion of an entry port on request</li> </ul>			
	Retrofitting (after technical clarification)	Retrofitting possible			
	Product requirements	Vacuum oven VD			
<b>Order no.</b>		N30m	N31m	N32m	N33m



## Entry port | Silicone

In addition to moisture tight entry seals, which are welded into the device, m360 offers plugged silicone entry ports in a moisture tight design.




Model sizes/Description			ø40	ø80	ø120
Description	Specification		Silicone entry ports for appliances with humidity control or special sealing requirements		
Technical data	Clear diameter	mm	40	80	120
	Technical feature		moisture-tight		
	Side panel		-		
	Back panel		●		
	Closure		Plug		
	Material		Silicone		
Ordering information	Further options (only after technical clarification)		Lateral mounting on request		
	Retrofitting (after technical clarification)		Retrofitting possible		
<b>Order no.</b>			N20m	N21m	N22m



## Entry ports

Round entry ports are mainly used for the insertion measuring sensors. The Memmert standard portfolio offers side entry ports with a diameter of up to 23 mm. m360 manufactures special entry ports with wider diameters that can also be individually positioned. Depending on the design, round entry ports can be closed with an ear cap or slider.




Model sizes/Description		Ø14	Ø23	Ø38	Ø57	Ø80	Ø100	Ø120						
Description	Specification	Entry ports for devices without special sealing requirements												
Technical data	Clear diameter	mm	14	23	38	57	80	100	120					
	Technical feature	Not welded = not moisture-tight												
	Side panel	-	●	-	●	-	●	-	●					
	Back panel	●	-	●	-	●	-	●	-					
	Closure	Flap				Slider								
Ordering information	Further options (only after technical clarification)	Intermediate sizes from Ø57mm possible												
	Additional information	Sealing plugs for entry ports available as accessories												
	Retrofitting (after technical clarification)	Retrofitting possible												
<b>Order no.</b>		D6m	N51m	F6m	F7m	N56m	F8m	N59m	N60m	N61m	F9m	N63m	N64m	N65m



### Modification | ESD protection


m360 offers an ESD appliance modification as an accessory. In addition to replacing the standard seal with an ESD seal, the ESD modification package offers further safety precautions against build-up of electrostatic charge in the appliance, such as ESD- capable castors and earthing bolts.



Model sizes/Description		110	160	260	450	750/1060
Description	Specification	ESD-safe equipment and connection options for earth connections				
	Design	Use of ESD-compatible components				
Technical data	Technical feature	ESD protection				
	Temperature resistance (additional option may be required)	°C	180			
Ordering information	Retrofitting (after technical clarification)	Retrofitting not possible				
	Restrictions	Appliance sizes without castors require an additional castor frame				
	Required option (depending on device type)				A8	R9
<b>Order no.</b>					M3m	

**Modification | Insulation protection**

The Memmert appliances are often used in cleanrooms. In order to reduce the particle load in the interior of the appliance, the rock wool insulation can also be packed in aluminium foil.




Model sizes/Description		Insulation protection
Description	Specification	Reduction of particle pollution through insulation wrapped in aluminium foil
Ordering information	Additional information	<ul style="list-style-type: none"> <li>No improvement in the ISO cleanroom class</li> <li>Heat output of the appliances deviates from the standard</li> </ul>
	Retrofitting (after technical clarification)	Retrofitting not possible
	Product requirements	Insulation material: rock wool/glass wool
<b>Order no.</b>		M2m



### Heavy-duty | Base frame


Heavy components can exceed the maximum load of the appliance due to their high weight. A heavy-duty base frame is therefore recommended for particularly heavy loads. Since the feet of the shelf unit stand directly on the foundation of the installation location, no load is placed on the appliance.



Model sizes/Description		260	450	750	1060
Description	Specification	Frame inside the interior, weight transfer to the floor.			
	Design	four entry ports on the base of the appliance			
Technical data	Position on the appliance	Interior			
	Technical feature	levelling feet			
	Max. Weight load frame	800 kg			
	Material	stainless steel 1.4301 (ASTM 304)			
Ordering information	Additional information	<ul style="list-style-type: none"> <li>Base frame for further heavy-duty interior fittings</li> </ul>			
	Retrofitting (after technical clarification)	Retrofitting not possible			
	Extensions	M11m M12m	M21m M23m	M22m M23m	M31m M32m
<b>Order no.</b>		M10m	M20m		M30m

### Heavy-duty | Grid (for rack)

Extension of the slide-in levels for the heavy-duty interior frame.




Model sizes/Description		260	450	750	1060
Description	Specification	Grid for heavy-duty interior frame			
Technical data	Position on the appliance	Interior			
	Max. Number of levels	number	3	6	
	Max. Weight load per level	kg	200		
	Material	stainless steel 1.4301 (ASTM 304)			
Ordering information	Retrofitting (after technical clarification)	Retrofitting not possible			
	Product requirements	Heavy-duty base frame with interior frame			
	Required option (depending on device type)	M10m M11m	M20m M21m	M20m M22m	M30m M31m
<b>Order no.</b>	M12m	M23m		M32m	



## Heavy-duty | Interior frame


The heavy-duty base frame of the appliance can be extended with a complete interior frame. Additional shelves can be used to increase the support surface in the interior considerably.



Model sizes/Description		260	450	750	1060
Description	Specification	Heavy-duty interior frame			
Technical data	Position on the appliance	Interior			
	Max. Number of levels	number	3	6	
	Max. Weight load per level	kg	200		
	Max. Weight load frame	kg	800		
	Material	stainless steel 1.4301 (ASTM 304)			
Ordering information	Additional information	incl. one heavy-duty grid			
	Retrofitting (after technical clarification)	Retrofitting not possible			
	Product requirements	Heavy-duty base frame			
	Required option (depending on device type)	M10m	M20m		M30m
	Extensions	M12m	M23m		M32m
<b>Order no.</b>		M11m	M21m	M22m	M31m

## Telescopic slides | Sizes 55 to 160

Telescopic slides enable Memmert appliances to be loaded more ergonomically. The pull-outs can be combined with grids and perforated shelves and each pull-out can carry a load of 60 kg. As a rule, three to four telescopic rails can be installed per appliance, but more pull-outs are possible after technical clarification.




Model sizes/Description		55/75	110/160	55/75	110/160		
Description	Specification	Telescopic slides					
	Design	Grid		Perforated sheet			
Technical data	Position on the appliance	Interior					
	Max. Weight load per level	60 kg					
	Pull-out length	A	mm	265	310	265	310
	Material	stainless steel 1.4301 (ASTM 304)					
Ordering information	Additional information	<ul style="list-style-type: none"> <li>Quantity and position only after technical clarification</li> <li>Mind the maximum weight load of the appliance</li> </ul>					
	Retrofitting (after technical clarification)	Retrofitting not possible					
	Restrictions	Levels cannot be modified or changed afterwards.					
<b>Order no.</b>		M6m		M7m			



### Door | Limit opening angle

This equipment enables efficient use of space thanks to customisable opening angles of 95° or 120°, allowing the working range to be more flexible and functional.




Model sizes/Description		all sizes	
Description	Specification	Limiting the door opening angle	
Technical data	Technical feature	95°	120°
Ordering information	Additional information	For appliances with double doors, the hinges can be combined as required	
	Retrofitting (after technical clarification)	Retrofitting possible	
<b>Order no.</b>		Q0m	Q1m



## Door | Replace door stop

For ergonomic reasons or to save space on-site, it may be advisable for some cases to mount the door hinge on the left-hand side so that the door handle is on the right-hand side. m360 attaches the door handle to your desired side at the factory.




Model sizes/Description		30 to 410
Description	Specification	Moving the door hinge to the other side (left - right)
Ordering information	Additional information	For appliances with a glass inner door, the stop is also replaced
	Retrofitting (after technical clarification)	Retrofitting possible
	Product requirements	single-leaf door
<b>Order no.</b>		B8m



## Door | Seals


For special applications, it may be necessary to replace the standard seal of a Memmert appliance with low-silicone seals. In this case, a low-silicone Viton seal is recommended, which is heat-resistant up to 250°C for long-term processes. ESD seals are conductive and therefore provide greater protection against electrostatic discharge.



Model sizes/Description		ESD	Viton	
Description	Specification	Replacing the standard seal		
	Design	ESD protection	Silicone-free Viton seal	
Technical data	Temperature resistance (additional option may be required)	°C	180	250
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible		
	Required option (depending on device type)	A8		
<b>Order no.</b>		Q3m	Q2m	

**Door | Handlebar**

Particularly for XXL appliances, but also for smaller Memmert appliances that are opened and closed very frequently, it is advisable to replace the aluminium handlebar with a reinforced handlebar made of stainless steel.




Model sizes/Description		from size 410
Description	Specification	Reinforced door handlebar
Technical data	Material	stainless steel 1.4301 (ASTM 304)
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible
	Product requirements	for appliances from size 410
<b>Order no.</b>		Q4m



## Door | Magnetic door lock

Certain processes require the appliance door to be frequently opened. For more ergonomic handling, a magnetic closure is the perfect solution. Thanks to the magnet, the door closes automatically and can withstand constant load of opening and closing.




### Model sizes/Description

Description	Specification	Door lock
	Design	magnetic
Technical data	Temperature resistance (additional option may be required) °C	200
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible
	Product requirements	single-leaf door
	Required option (depending on device type)	A8
<b>Order no.</b>		M0m

### Door | Full-sight glass

The full-sight glass door provides a direct view of the interior without having to open the oven. It is the ideal solution for laboratories and industry, offering maximum control over work processes.




Model sizes/Description		110	260	410	750/1060	410	750/1060
Description	Specification	Outer door inspection window					
	Design	Unheated				Heated	
Technical data	Inspection window width	mm	388		241		
	Inspection window height	mm	425	744	1144		
Ordering information	Additional information	Inner glass door not applicable					
	Retrofitting (after technical clarification)	Retrofitting possible					
<b>Order no.</b>			B0m			B3m	



## Maintenance | Electronics box

To save space or for certain applications in cleanrooms, appliances are often recessed into a wall. As a result, the cover of the electronics box can no longer be opened. With the „maintenance hatch“ modification, the conventional appliance cover is replaced by a two-part cover.




### Model sizes/Description


Description	Specification	Maintenance hatch		
Technical data	Position on the appliance	Electronic box		
	Distance front panel	mm	237	
	Height flap	mm	350	450
	Material	textured stainless steel 1.4016 (ASTM 430)		
Ordering information	Retrofitting (after technical clarification)	Retrofitting not possible		
	Product requirements	for appliances up to size 410		
<b>Order no.</b>		M1m		

### Supply air | Filter system

Supply air filter system for controlled air quality in the interior. Depending on the version, the air filter holders offer capacity for up to two filters. Depending on the requirements, the filter holder can be fitted with a suitable filter.



Model sizes/Description		Two Filter system	Single Filtersystem
Description	Specification	Filter holder for supply air opening	
	Design	2 Filter	1 Filter
Technical data	Material	stainless steel 1.4301 (ASTM 304)	
Ordering information	Additional information	Filter holder without filter	
	Retrofitting (after technical clarification)	Retrofitting possible	
	Required option (depending on device type)	R9	
	Extensions	E07439 E07438	
<b>Order no.</b>		R82m	R81m



**TECHNICAL  
INFORMATIONS  
ELECTRONICS**




Exhaust air flap | Rear panel  
Compressed air drying  
External control cabinet  
External control unit  
Light | Lighting interior  
Light | Lighting interior  
Light | Lighting LED light strips - grid  
Light | Lighting LED light strips - K1 extension  
Exhaust air extraction | Rapid cooling fan  
Sensors | Temperature  
Sensors | Further parameters  
Interfaces | Potential free contact  
Interfaces | 4-20mA  
Socket outlet | Interior  
Signal lamp | Top cover  
Signal lamp | Control panel  
Door | Safety interlock  
Door | CO<sub>2</sub> incubator locking



## Exhaust air flap | Rear panel


The exhaust air flap enables precise regulation of the air exchange, resulting in an increased air exchange rate to ensure optimal air quality.



Model sizes/Description		Exhaust air flap
Description	Specification	Electronically adjustable exhaust air flap
Technical data	Technical feature	Opening degree adjustable in 10% increments
	Option to control or programme the parameters	AtmoCONTROL AtmoWEB ControlCOCKPIT
Ordering information	Additional information	Programmable via AtmoCONTROL or AtmoWEB
	Retrofitting (after technical clarification)	Retrofitting possible
<b>Order no.</b>		E6m

## Compressed air drying

This option extends the temperature-humidity operational range of the appliance. Thanks to the efficient dehumidification of the interior by means of a compressed air supply, it is possible to dehumidify the air in the appliance to a relative humidity of up to 10 % even at 10 °C.




Model sizes/Description		Compressed air drying	
Description	Specification	Efficient dehumidification of the interior using a compressed air supply (e.g. +10 °C with 10% rh)	
Technical data	Working pressure depending on operating condition	bar	0.2
	Max. permissible pressure	bar	5
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible	
	Required option (depending on device type)	U1m	
<b>Order no.</b>		C9m	



## External control cabinet


This option offers a convenient and device-independent arrangement of all electronic components and the control unit, so that operation and maintenance can be carried out from a different location.



Model sizes/Description		External control cabinet
Description	Specification	External units for electronic components
Technical data	Protection class	IP20
	Technical feature	Connection cable 3m
Ordering information	Further options (only after technical clarification)	Alternative lengths of connection cables possible
	Retrofitting (after technical clarification)	Retrofitting not possible
<b>Order no.</b>		E5m

## External control unit

Equipping the Memmert oven with an external control unit allows the oven to be conveniently controlled remotely from any other position in the room. This ensures flexible and efficient operation.




Model sizes/Description		External control unit
Description	Specification	External unit for regulator module
Technical data	Protection class	IP20
	Technical feature	Connection cable 3m
Ordering information	Further options (only after technical clarification)	Alternative lengths of connection cables possible
	Retrofitting (after technical clarification)	Retrofitting possible
<b>Order no.</b>		E4m



## Light | Lighting interior


To provide a better view of the load inside the appliance chamber, m360 offers special lighting for monitoring the load. The lights can be switched on manually as well as being switched on continuously.



Model sizes/Description		110 - 260	450 - 750
Description	Specification	Interior lighting for observation	
	Design	1 light	2 lights
Technical data	Position on the appliance	Interior	
	Technical feature	switchable: ON/OFF	
	Option to control or programme the parameters	AtmoCONTROL AtmoWEB ControlCOCKPIT	
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible	
<b>Order no.</b>		R0m	

## Light | Lighting interior

Thanks to the special m360 solution, lighting is now available for the vacuum oven. The light replaces the thermoshelf and illuminates the interior of the vacuum chamber to ensure optimum viewing of the process from the outside.




Model sizes/Description		Vacuum oven VO
Description	Specification	Interior lighting for observation
Technical data	Position on the appliance	Flange socket
	Option to control or programme the parameters	AtmoCONTROL AtmoWEB ControlCOCKPIT
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible
	Product requirements	Vacuum oven VO
	Restrictions	at least one slot for a thermoshelf is not available (depending on product size)
<b>Order no.</b>		ROm



## Light | Lighting LED light strips - grid

For certain applications, such as plant or mushroom cultivation, it is vital to illuminate the load from above. The m360 LED grid offers this feature. The light can be dimmed in 1% increments and is available in warm white, cool white and a combination of both. Certain sequences can be programmed via AtmoCONTROL, e.g. a day/night simulation.




Model sizes/Description		Cold white + warm white		Cold white		Warm white		
Description	Specification	Dimmable LED light strips (0-100% in 1% increments)						
	Design	Mounting on the grid						
Technical data	Position on the appliance	Grid	Reinforced grid	Grid	Reinforced grid	Grid	Reinforced grid	
	Colour temperature	6500 + 2700		6500		2700		
	Protection class	IP44						
	Limited working temperature range	°C +15 to +40						
	Option to control or programme the parameters	AtmoCONTROL AtmoWEB ControlCOCKPIT						
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible						
	Product requirements	TwinDISPLAY						
<b>Order no.</b>		J18m	J28m	J17m	J27m	J19m	J29m	
LED equipment per model size		110	260	410	750	1060	1400	2200

Equipment	Max. Number of illuminated levels	number	2	3		5	6
	LED light strips per slide-in level	number	5		3		



## Light | Lighting LED light strips - K1 extension

The m360 solution combines LED lighting with the storage of heavy loads. Mounted LED light strips on the K1 heavy-duty extension provide light-intensive illumination and enable precise control in 1% increments. Available as warm white, cool white or combined. Programmable sequences with AtmoCONTROL, such as day/night simulation.




Model sizes/Description		Cold white	Warm white	Cold white + warm white
Description	Specification	Dimmable LED light strips (0-100% in 1% increments)		
	Design	Mounting on the K1 extension		
Technical data	Colour temperature	6500	2700	2700 + 6500
	Protection class	IP44		
	Limited working temperature range °C	+15 to +40		
	Option to control or programme the parameters	AtmoCONTROL AtmoWEB ControlCOCKPIT		
Ordering information	Additional information	K1 extension only required for appliance sizes 410 and 750		
	Retrofitting (after technical clarification)	Retrofitting possible		
	Product requirements	TwinDISPLAY		
	Required option (depending on device type)	K1		
<b>Order no.</b>		J07m	J09m	J08m

LED equipment per model size		410	750	1060	1400	2200
Equipment	Max. Number of illuminated levels	number		7	8	9
	LED light strips per slide-in level	number		2		



## Exhaust air extraction | Rapid cooling fan


A rapid cooling fan has two functions: On the one hand, it cools the appliance many times faster after a tempering process without having to open the door. On the other hand, the fan increases the air exchange rate of the appliance. The rapid cooling fan can be adjusted in 10 per cent increments depending on the exhaust air flap and can be controlled via a controller on the exhaust air flap, but also via the AtmoCONTROL software.



Model sizes/Description		Rapid cooling fan
Description	Specification	Efficient extraction of exhaust air
	Design	Adjustable via exhaust air flap in 10% increments
Technical data	Position on the appliance	rear position
	Technical feature	Hose connection $\varnothing 60.3\text{mm}$
	Option to control or programme the parameters	AtmoCONTROL AtmoWEB ControlCOCKPIT
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible
	Product requirements	Exhaust air opening
	Restrictions	Performance data deviate from the standard appliances
<b>Order no.</b>		E3m

## Sensors | Temperature

With additional, freely positionable temperature sensors, users have the option of measuring not only the air temperature, but also the temperature of the charging material or water in a tray, for example.




Model sizes/Description		Flexible temperature sensors			
Description	Specification	PT100 can be flexibly positioned in the interior			
	Design	Connection socket (4-pin) on the rear of the appliance	Temperature indication on the display		
Technical data	Max. Quantity	number	1	2	3
	Option to evaluate the parameters		Interface for external systems	AtmoCONTROL AtmoWEB ControlCOCKPIT	
Ordering information	Retrofitting (after technical clarification)		Retrofitting possible		
<b>Order no.</b>		H4m	H8m	H8m2	H8m3



## Sensors | Further parameters


In special applications, it is necessary to monitor a parameter such as relative humidity, although the device cannot actively control the parameter. In this case, m360 can, for example, install a humidity sensor in a UF device, which displays and logs the relative humidity. However, it is still not possible to control the humidity parameter. The m360 team of experts will be happy to advise you.



Model sizes/Description		CO2	Humidity	O2
Description	Specification	Measuring sensors for non-adjustable parameters		
Technical data	Position on the appliance	Interior		
	Technical feature	CO <sub>2</sub>	Feuchte	O <sub>2</sub>
	Measuring range of the sensor	0 to 20 %CO <sub>2</sub>	0 to 100 %rh	0 to 25 %O <sub>2</sub>
	Option to evaluate the parameters	AtmoCONTROL AtmoWEB ControlCOCKPIT		
Temperature resistance (additional option may be required)	°C	70	180	200
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible		
	Product requirements	TwinDISPLAY		
	Restrictions	Parameters are displayed and logged on the appliance display. They cannot be controlled.		
	Required option (depending on device type)	A8		
<b>Order no.</b>		U2m	U1m	U0m

## Interfaces | Potential free contact

Potential-free contacts are available for all UF/IF appliances in the Memmert standard portfolio. The m360 team of experts has developed a special configuration for the other product models. Using the AtmoCONTROL appliance software, a dependency can be programmed for options H7 / H72/ H74 and an additional signal is sent to the potential-free contact, e.g. the command to activate ventilation when 300 °C is reached. With the H5option, an interface is mounted on the rear panel of the appliance where external devices, e.g. for humidity monitoring, can be connected.




Model sizes/Description		Potential free contact					
Description	Specification	Potential free switching contact (24V / 2A)					
	Design	programmable				Target value reached	
Technical data	Max. Quantity	number	1	2	3	4	1
	Option to control or programme the parameters		AtmoCONTROL				-
	Option to evaluate the parameters		AtmoWEB Interface for external systems				Interface for external systems
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible					
	Product requirements		TwinDISPLAY				-
<b>Order no.</b>		H7m	H7m2	H7m3	H7m4	H5m	



## Interfaces | 4-20mA


To log the relative humidity in the appliance via an external system, a 4-20mA interface can be installed to output the actual value.



Model sizes/Description		4-20mA interface
Description	Specification	Parameter evaluation humidity via 4-20mA
	Design	0 to 100% rh = 4-20mA
Technical data	Technical feature	Humidity
	Option to evaluate the parameters	Interface for external systems
Ordering information	Additional information	An additional option may be required for devices without a humidity sensor
	Retrofitting (after technical clarification)	Retrofitting possible
<b>Order no.</b>		V7m

## Socket outlet | Interior

Many users connect additional appliances to their Memmert oven, e.g. shaking devices or scales. If you cannot have a entry port installed and would like more than one interior socket for your appliance, you may contact m360.




Model sizes/Description		Socket outlet			
Description	Specification	Socket outlet in the interior			
	Design	switched on/off via toggle switch on the control panel		switched off via the main switch	
Technical data	Position on the appliance	Interior			
	Socket outlet for plug type	Type F „Schuko-plug“	Type G „UK-plug“	Type F „Schuko-plug“	Type G „UK-plug“
	Protection class	IP68	IP66	IP68	IP66
	Technical feature	230 V / 2,2 A			
	Temperature resistance (additional option may be required)	°C	60	40	60
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible			
	Required option (depending on device type)	A8			
<b>Order no.</b>		R4m	R44m	R3m	R33m



## Signal lamp | Top cover

Signal lights allow the user to recognise the ongoing process in the Memmert appliance: green for operation, yellow for standby and red for an error message. They are particularly useful when several appliances are in use at the same time. They are mounted on the upper panel of the appliance.




Model sizes/Description		Signal lamp   Top cover	
Description	Specification	Status signal lamps <ul style="list-style-type: none"> <li>• Red = error/alarm</li> <li>• Yellow = standby</li> <li>• Green = programme is running</li> </ul>	
Technical data	Position on the appliance	Top cover	
	Technical feature	single-colour lights	multicolor light
	Max. Quantity	number	3
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible	
	Product requirements	TwinDISPLAY	
	Restrictions	Other options may be excluded by the installation	
	Required option (depending on device type)	H6 H72	
<b>Order no.</b>		E1m	E2m



## Signal lamp | Control panel

Signal lights allow the user to identify the current status of the ongoing process in the Memmert appliance: green for ongoing operation, yellow for standby and red for an error message. They are particularly useful when several appliances are in use at the same time. They are mounted on the upper panel of the appliance.




Model sizes/Description		Signal lamp   Control panel
Description	Specification	Status signal lamps <ul style="list-style-type: none"> <li>• Red = error/alarm</li> <li>• Yellow = standby</li> <li>• Green = programme is running</li> </ul>
	Technical data	Front cover
Technical data	Position on the appliance	single-colour lights
	Technical feature	
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible
	Product requirements	TwinDISPLAY
	Restrictions	Other options may be excluded by the installation
	Required option (depending on device type)	H6 H72
<b>Order no.</b>		E0m



## Door | Safety interlock


The safety door interlock meets high safety standards and can be controlled via external systems. Reliable safety switches make it possible to control almost all process steps and ensure process safety inside the appliance.




Model sizes/Description		Safety door interlock
Description	Specification	Electronic door lock with safety switch: <ul style="list-style-type: none"> <li>• Safety switch with guard locking and integrated evaluation electronics                             <ul style="list-style-type: none"> <li>• Locking force up to 6500 N</li> <li>• Short circuit monitoring</li> </ul> </li> <li>• 2 safety outputs (semiconductor outputs)</li> <li>• Category 4 / PL e according to EN ISO 13849-1                             <ul style="list-style-type: none"> <li>• With M23 plug connector (RC18)</li> </ul> </li> </ul>
	Technical data	
Ordering information	Option to control or programme the parameters	Interface for external systems
	Option to evaluate the parameters	Interface for external systems
	Additional information	Control and programming is performed by external systems
	Retrofitting (after technical clarification)	Retrofitting possible
	Product requirements	single-leaf door
<b>Order no.</b>		E7m

## Door | CO<sub>2</sub> incubator locking

An addition to the standard option portfolio for the ICO appliance series is now available. The process-dependent door lock, developed by m360, can now also be used for CO<sub>2</sub> incubators. Process security and access control can be programmed via AtmoCONTROL or directly on the ControlCOCKPIT.



Model sizes/Description		CO <sub>2</sub> incubator door locking
Description	Specification	Electronic door lock
	Design	Process-dependent programmable door lock
Technical data	Option to control or programme the parameters	AtmoCONTROL AtmoWEB ControlCOCKPIT
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible
	Product requirements	Incubator ICO
<b>Order no.</b>		D4m



**TECHNICAL  
INFORMATIONS  
ACCESSORY**

Setup | Stacking rack for vacuum ovens

Setup | Stacking frame

Plug | Stopper entry port

Exhaust air | Condensate drainage

Connection cable | lengths

Shelves | Perforated sheets, grids


Telescopic slides | Sizes 260 to 1060

Heavy-duty | Floor grid

Supply air | Filter

## Setup | Stacking rack for vacuum ovens


The stacking rack for vacuum ovens enables the space-saving and ergonomic installation of various sizes of vacuum ovens. This allows work processes to be optimised and organisation in the laboratory to be flexible.



Model sizes/Description		29/49	101
Description	Specification	Rack for stacking vacuum ovens VO.	
	Design	Stacking of sizes 29 and 49	Stacking of size 101
Ordering information	Additional information	<ul style="list-style-type: none"> <li>For appliances of different sizes, the larger appliance must always be at the bottom.                             <ul style="list-style-type: none"> <li>Stacking of different sizes possible.</li> </ul> </li> </ul>	
	Retrofitting (after technical clarification)	Retrofitting possible	
<b>Order no.</b>		B50577	B59029

## Setup | Stacking frame

To save space in the laboratory or workshop, Memmert appliances can be stacked. The m360 stacking frame enables appliances of different sizes to be stacked safely, e.g. a UF55 on top of an IF160. For safety reasons, only two appliances can be stacked on top of each other, with the larger of the two appliances at the bottom.




Model sizes/Description		IPPeco /HPPeco			U / I	
Description	Specification	Stacking frame for stacking appliances.				
	Design	110	260	75	160	260
Technical data	Material	stainless steel 1.4301 (ASTM 304)				
Ordering information	Additional information	<ul style="list-style-type: none"> <li>For appliances of different sizes, the larger appliance must always be at the bottom.                             <ul style="list-style-type: none"> <li>Stacking of different sizes possible.</li> </ul> </li> </ul>				
	Retrofitting (after technical clarification)	Retrofitting possible				
<b>Order no.</b>		B57998	B58000	B50598	B50599	B50600



## Plug | Stopper entry port

In order to keep the climatic conditions of an appliance with a entry port constant during operation, the use of a stopper is recommended.



Model sizes/Description				Ø14	Ø23	Ø37	Ø57	Ø80	Ø100	Ø120
Description	Specification	Stopper for entry ports								
Technical data	Outer diameter	A	mm	13	23	38	55	75	96	116
	Outer diameter	B	mm	17.5	30	42	66	87	108	128
	Length	C	mm	25	33	46	52			
	Temperature resistance (additional option may be required)		°C	180						
Material				Silicone (sponge)						
<b>Order no.</b>				E05790	E05361	E05362	E07865	E05625	E07866	E05920


Model sizes/Description				Ø23	Ø37	Ø57	Ø100
Description	Specification	Stopper for entry ports					
Technical data	Outer diameter	A	mm	21	37	50.5	94
	Outer diameter	B	mm	27	45	59.5	107
	Length	C	mm	30	42	45	64
	Temperature resistance (additional option may be required)		°C	200			
Material				Silicone			
<b>Order no.</b>				E05387	E07479	E05389	E05386

Model sizes/Description				Ø23	Ø37	Ø57	Ø80	Ø100-120
Description	Specification	Stopper for entry ports						
Technical data	Outer diameter	A	mm	21	33.2	53.8	69.9	88.9
	Outer diameter	B	mm	27	41.2	63.7	88.9	127
	Length	C	mm	28.3	25.4		38.1	50.8
	Temperature resistance (additional option may be required)		°C	300				
Material				Silicone				
<b>Order no.</b>				E07698	E07699	E07700	E07701	E07702



## Exhaust air | Condensate drainage


The drying process in a Memmert appliance can produce humid air that needs to be removed. The moist air may condense in the exhaust air and drip down. A special valve with hose nozzle allows the water to run off in a controlled manner. The exhaust air condensate trap is made of stainless steel and is very easy to install - even retrospectively and without a technician.



Model sizes/Description		Condensate drainage
Description	Specification	Condensate drainage for exhaust air opening <ul style="list-style-type: none"> <li>Exhaust air pipe <math>\varnothing 60,3</math> mm</li> <li>Condensate drainage <math>\varnothing 18</math> mm</li> </ul>
Technical data	Material	stainless steel 1.4301 (ASTM 304)
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible
<b>Order no.</b>		B37684

## Connection cable | lengths

With a generous length of 5 metres, the power connection cable for Memmert ovens (230V) ensures a convenient connection and enables the appliance to be supplied with power if it cannot be set up in the immediate vicinity of a socket outlet.




### Model sizes/Description

Description	Specification	Extended power connection cable
	Design	5 m
Technical data	Socket outlet for plug type	Type F „Schuko-plug“
	Protection class	IP20
	Technical feature	230 V
<b>Order no.</b>		B58607

## Shelves | Perforated sheets, grids

The reinforced perforated shelves are equipped with additional reinforcements and are designed for higher loads. The load capacity increases to 80 kg for appliances without interior fittings (K1) and to 100 kg for appliances with interior fittings (K1). It is essential to ensure that the maximum load limit of the appliance is not exceeded.




Model sizes/Description		260 / 410	450 / 750	1060	260 / 410	450 / 750	1400 / 2200	
Description	Specification	Reinforced shelves and special sizes						
	Design	Perforated sheet					Reinforced grid	
Technical data	Max. Weight load per level	kg	100		80	30	60	
	Material	stainless steel 1.4301 (ASTM 304)						
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible						
	Restrictions	Mind the maximum weight load of the appliances.						
	Required option (depending on device type)		K1			-		
<b>Order no.</b>		B58460	B58461	B58462	B57086	B52092	B59007	B59004



## Telescopic slides | Sizes 260 to 1060

Telescopic slides enable a Memmert appliance to be loaded more ergonomically. The pull-outs can be combined with grids and perforated shelves and each pull-out can carry a load of 60 kg. As a rule, three to four telescopic slides can be installed per appliance, but more slides are also possible after technical clarification.




Model sizes/Description		260/410		450/750		1060	
Description	Specification	Installation of telescopic slides					
	Design	Grid	Perforated sheet	Grid	Perforated sheet	Grid	Perforated sheet
Technical data	Position on the appliance	Interior K1 extension					
	Max. Weight load per level	kg		60			
	Pull-out length	A	mm	410	505	660	
	Material	stainless steel 1.4301 (ASTM 304)					
Ordering information	Additional information	Number of telescopic slides after clarification					
	Retrofitting (after technical clarification)	Retrofitting not possible					
	Required option (depending on device type)	K1					
<b>Order no.</b>	B43102		B51618	B43093	B44476	B43100	B56888

### LED equipment per model size

Technical data	Max. Levels with telescopic runners	number	5	7	4	7
----------------	-------------------------------------	--------	---	---	---	---

### Heavy-duty | Floor grid

In some applications, the load exceeds the permissible maximum weight of 25 kg - 30 kg per shelf level. The usual perforated shelves and grids must then no longer be loaded. In such cases, a heavy-duty floor grid is recommended that can carry loads of up to 300 kg, depending on the appliance. The heavy-duty floor grid is made of stainless steel and can be combined with all options. On request, the floor grating can also be adapted so that special interior fittings can be attached.




Model sizes/Description		30	55/75	110/160	256	260/410	450/750	1060
Description	Specification	Floor grid for product storage at the base of the appliance						
Technical data	Position on the appliance	Below the appliance						
	Max. Number of levels	number	1					
	Material	stainless steel 1.4301 (ASTM 304)						
Ordering information	Restrictions	The maximum weight load of the appliance must not be exceeded						
<b>Order no.</b>		E50944	E46407	E36692	E41309	E36704	B31972	E40438



### Supply air | Filter

G3 filter (option R8) with 80% separation efficiency can be configured for Memmert appliances from series production. This filter prevents coarse particles from entering the interior of the appliance. For users who want to work with as few particles as possible but do not have a clean room, m360 offers the HEPA filter H13 (filtration efficiency of 99.95%). Other filter types, such as an F7 ultra-fine filter, are also available on request and after technical consultation.



#### Model sizes/Description

Description	Specification	Filter for supply air filter holder	
Technical data	Filter class	F7	H13
Ordering information	Retrofitting (after technical clarification)	Retrofitting possible	
	Required option (depending on device type)	R81m R82m R8	
<b>Order no.</b>		E07438	E07439

**PERSONAL NOTES**

A large rectangular area filled with horizontal dotted lines, intended for taking personal notes. The lines are evenly spaced and extend across the width of the page.

# SOLUTIONS FOR YOU



Memmert GmbH + Co. KG  
Äußere Rittersbacher Straße 38 | 91126 Schwabach  
Tel. +49 9122 925 -0 | Fax +49 9122 14585  
E-Mail: [m360@memmert.com](mailto:m360@memmert.com)

[www.memmert.com/m360](http://www.memmert.com/m360)  
[facebook.com/memmert.family](https://facebook.com/memmert.family)  
[instagram.com/memmert.family](https://instagram.com/memmert.family)  
[linkedin.com/company/memmert-gmbh-co-kg](https://linkedin.com/company/memmert-gmbh-co-kg)